

COMFORT GLIDE



Product Instructions

A Comfort Glide Repositioning Sheet **B** Comfort Glide Wedges **C** Comfort Glide Drypad

How to use the Comfort Glide Sheet and Drypad



1 Unfold and position the glide sheet next to the patient. Be sure to align the top of the sheet with the patient's shoulders.



2 Place a Comfort Glide Drypad over the glide sheet, printed side down.



3 Fold the edges of both lengthwise to tuck under the patient.



4 Roll the patient in the opposite direction.



5 Unfold the glide sheet and pad. Gently return the patient to his/her back.



6 Using the glide sheet handles, gently slide the patient as needed to centre him/her on the bed.

How to use Comfort Glide Wedges



1 Gently roll the patient onto his/her side.



2 Locate the patient's sacrum. Place one wedge two inches above the sacrum.



3 Place another wedge two inches below the sacrum.



4 Make sure the wedges are about one fist apart.



5 Slowly roll the patient onto the wedges. Be sure the sacrum is 'floating' between the two wedges.

Caution:

1. Comfort Glide components are not designed for patient lifting.
2. These products are not to be used in combination with other lifting or positioning devices.
3. For best outcomes, only use a combination of Comfort Glide Sheets, Wedges and Drypads.
4. ATTENTION: If using wedges with Comfort Glide, do not use a Comfort Glide Sheet to slide the patient onto the wedges. Use the log-roll method to place the wedges under the patient.
5. Always refer to and follow your facility's guidelines for proper patient handling.

For additional information on Comfort Glide, please contact your account manager or visit us at: www.medline.eu/uk