


CLOSED WOUND DRAINAGE SYSTEMS

- en** INSTRUCTIONS FOR USE (p.2-3)
- de** GEBRAUCHSANWEISUNG (p.4-5)
- fr** NOTICE D'UTILISATION (p.6-7)
- nl** GEBRUIKSAANWIJZING (p.8-9)
- it** ISTRUZIONI PER L'USO (p.10-11)
- es** INSTRUCCIONES DE USO (p.12-13)
- pt** INSTRUÇÕES PARA USO (p.14-15)
- sv** INSTRUKTIONER FÖR ANVÄNDNING (p.16-17)
- fi** KÄYTTÖOHJEET (p.18-19)
- da** INSTRUKTION (p.20-21)
- no** BRUKSANVISNING (p.22-23)

		
REF DYNJWE1305 (100 cc/mL)	REF DYNJWE2000 (200 cc/mL)	REF DYNJWE1000 (400 cc/mL)

en CLOSED WOUND DRAINAGE SYSTEMS

Description

Medline closed wound drainage systems are single use sterile products. The closed wound drainage system consists of an internal drain and a bulb shaped fluid collection reservoir. The drain is either perforated, channeled, or ducted in nature and is manufactured to be radio-opaque. The wound drain is surgically placed into the surgical wound site using a trocar. Exudate is collected through a wound drain into a fluid collection Reservoir.

Indications

The drains are commonly used post-operatively for collecting bodily fluids from surgical sites.

Warnings/Precautions

1. Drains, Trocars, and Reservoirs are for single use only and should be discarded after use.
2. Drains or drain tubing should not be handled with any instruments. Sharp objects can lead to tearing or weakening and breakage of the drain.
3. Do not use if package has been opened or if the seal has been compromised.
4. In order to allow for easy removal; the drain and tubing should not be curled, stretched or sutured into place. Drains should be placed and removed with care avoiding fast movements applying slow and steady pressure.
5. Utilize sterile technique when making drain connections.
6. The product can be used with medical grade adhesive tape.
7. During drain placement, take care to not cut, tear or damage the drain as this can lead to tearing or complete breakage.
8. Careful attention should be directed to make sure the drains do not clog or become occluded when they are still in use. If occlusion occurs, suction directly to the drain itself can be administered to relieve occlusion. The reservoir bulb should not be allowed to fill completely, this can lead to increased likelihood of back-flow contamination.
9. All drainage holes, channels, or ducts must be located within the wound site creating an air tight seal to maintain vacuum. An air tight seal must be maintained to prevent the reservoir from filling with air which can increase the likelihood of backflow.
10. Do not resterilize the product.

Contraindications

1. Do not suture through, cut, nick, tear or damage the drain as this may lead to breakage.
2. Do not use with any other suction vacuum source other than the 100 cc/mL, 200cc/mL or 400 cc/mL reservoir.

Risks/Side Effects

Possible risks include infection and inflammation. Prior to the surgery, the patient should be informed of such risks and complications that may occur depending on the patients degree of intolerance to any object placed into the body. A drain left in a patient for any period of time allows for tissue ingrowth around and into the drain which can cause breakage and discomfort upon removal.

Reservoir Initial Testing

For Bulb Reservoir (Medline), please follow the instructions / diagram below. For other vendor reservoir, please refer to the vendor's instructions.

1. Remove Reservoir (D) from packaging and test for efficacy as follows:
 - a. Fully compress Reservoir (D).
 - b. Close Disposal Port (C and A).
 - c. Occlude Patient Port (B and E) and release Reservoir (D).
 - d. If vacuum is effective, Reservoir (D) should remain compressed. If Reservoir (D) expands, it should be discarded.


Drain Connection/Drainage Placement

Place drain in wound site, ensuring that all holes, channels, or ducts are located within the wound site.


1. Attach the drain or adaptor to the reservoir via Patient Port (B).
2. For Dual drain connections (400 cc/mL bulb only): attach second drain or adaptor to Patient Port 2 (E).
3. Firmly compress Reservoir (D) to evacuate air.
4. Close Disposal Port (C) with Screw Cap (A).
5. Release Reservoir (D) – the closed wound drainage system is now active.

Reservoir should be located at patient level or below drain exit for adequate suction to be maintained. Drainage should be monitored such that Reservoir is not filled to excess.


Silicone Reservoir 100 cc/mL

Silicone Reservoir 200 cc/mL

Silicone Reservoir 400 cc/mL

Disposal

Please refer to your hospital policy for correct disposal of this device and its contents.

Risks Associated with Reuse

The wound drainage system devices are not validated for reuse. The design of the device does not facilitate cleaning and can result in an increased risk of infection if not properly decontaminated. The device may not perform as intended by the manufacturer if reused, resterilized, or decontaminated. Reuse of the product may result in material degradation resulting in a device failure.

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Single use only

Product Code

Manufactured by:  **Conod Medical Co., Limited**, No. 38 Hongfeng Road,
Baimao, Industrial Park, Guli Town, Changshu City,
Jiangsu Province, 215532. P.R. China

CE 0123 **EC REP** Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg Germany

Made in China- Hergestellt in China - Fabriqué en Chine - Geproduceerd in China - Fabricato in Cina -
Fabricado en China - Fabricado na China - Tillverkad i Kina - Valmistettu Kiinassa - Fremstillet i Kina - Laget i Kina.

Distributed in USA by: Medline Industries, Inc., Three Lakes Dr., Northfield, IL 60093 USA.

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

Distributed in Europe by: Medline International France, 5 Rue Charles Lindbergh, 44110, Châteaubriant, France

Medline is a registered trademark of Medline Industries, Inc.

RE17CSI 2017-05