en

## **GEMINI STERILIZATION WRAPS**

INSTRUCTIONS FOR USE (p.2-13)

de

## **GEMINI STERILISATIONSVLIESE**

GEBRAUCHSANWEISUNG (P.14-26)

fr

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nl

## **GEMINI STERILIZATIEVERPAKKINGEN**

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it

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es

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INSTRUCCIONES DE USO (P.65-77)

pt

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sv

## **GEMINI PACKSKYNKE**

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fi

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da

## **GEMINI STERILEMBALLAGE**

INSTRUKTION (P.115-126)

no

## **GEMINI STERILISERINGSOMSLAG**

BRUKSANVISNING (P.127-138)

cs

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ro

## **PUNGI DE STERILIZARE GEMINI**

INSTRUCTIUNI DE FOLOSIRE (P.151-164)

hu

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MODELS - MODELLE - MODÈLES - MODELLEN - MODELLI - MODELOS - MODELOS -  
MODELLER - MALLIT - MODELLER - MODELLER - MODELY - MODELE - MODELLEK:

GEM01XX-EU	-	GEM01XXT-EU	-
GEM11XX-EU	GEM11XXINT-EU	GEM11XXT-EU	GEM11XXS-EU
GEM21XX-EU	GEM21XXINT-EU	GEM21XXT-EU	GEM21XXS-EU
GEM31XX-EU	GEM31XXINT-EU	GEM31XXT-EU	GEM31XXS-EU
GEM41XX-EU	GEM41XXINT-EU	GEM41XXT-EU	GEM41XXS-EU
GEM51XX-EU	-	GEM51XXT-EU	GEM51XXS-EU

## PRODUCT DESCRIPTION

Gemini Sterilization Wrap is offered to the marketplace as bulk packages of single sheets, where in accordance with standard hospital practices, two sheets are then used to wrap a medical device or a collection of medical devices for sterilization. Gemini Bonded Sterilization Wrap comprises of two sheets of Gemini Sterilization Wrap ultrasonically sealed on two sides. This allows for convenient wrapping with two sheets simultaneously. Gemini Interleaved sheets bring both our multi-layer technology and a colour differentiation to sterilization units where sequential wrapping is in place.

The sheets of sterilization wrap are square or rectangular sheets of fabric produced using a five layer SSMMS (spunbond-spunbond-meltblown-meltblown-spunbond) process. The wrap is available in single, bonded or interleaved sheet configuration. The single wrap is blue. The bonded wrap is single colour (blue sheet bonded with a blue sheet) or two colour (blue sheet bonded with a pink sheet). The interleaved wrap is comprised of two single sheets of blue and pink. The wrap allows for aseptic opening of the sterilized package.

Gemini Sterilization Wraps are available in sizes ranging from 46x46 cm to 137x183 cm across six different material weights/models listed in table 1.

**Table 1. Gemini Sterilization Wrap Weights/Models**

Weight	Single Ply	Interleaved	Bonded Two Colour	Bonded Single Colour
Super Light Weight	GEM01XX-EU	Not applicable	GEM01XXT-EU	Not applicable
Light Weight	GEM11XX-EU	GEM11XXINT-EU	GEM11XXT-EU	GEM11XXS-EU
Regular Weight	GEM21XX-EU	GEM21XXINT-EU	GEM21XXT-EU	GEM21XXS-EU
Medium Weight	GEM31XX-EU	GEM31XXINT-EU	GEM31XXT-EU	GEM31XXS-EU
Heavy Weight	GEM41XX-EU	GEM41XXINT-EU	GEM41XXT-EU	GEM41XXS-EU
Super Heavy Weight	GEM51XX-EU	Not applicable	GEM51XXT-EU	GEM51XXS-EU

# INDICATIONS FOR USE

Gemini Sterilization Wraps are intended to be used to enclose another medical device that is to be sterilized by a healthcare provider using the following sterilization methods.  
The dry times below are “minimum” dry times, and will need to be adjusted in the real environment.

**Table 2. Pre-vacuum steam at 134°C from 3 to 18 minutes with the following dry times by weight:**

GEM0 – Super Light Weight	20 Minutes
GEM1 – Light Weight	20 Minutes
GEM2 – Regular Weight	20 Minutes
GEM3 – Medium Weight	20 Minutes
GEM4 – Heavy Weight	30 Minutes
GEM5 – Super Heavy Weight	30 Minutes

**Table 3. Gravity steam at 121°C from 15 to 30 minutes with the following dry times by weight:**

GEM0 – Super Light Weight	20 Minutes
GEM1 – Light Weight	20 Minutes
GEM2 – Regular Weight	20 Minutes
GEM3 – Medium Weight	20 Minutes
GEM4 – Heavy Weight	30 Minutes
GEM5 – Super Heavy Weight	30 Minutes

**Table 4. Ethylene Oxide (ETO Sterilization)**

60 Minutes exposures to 725mg/L of 100% ETO at 55°C and 40%-80% relative humidity followed by 12 hours aeration time.	Not applicable for GEM0.
---	--------------------------

**Table 5. Advanced Sterilization Products STERRAD® Sterilization System – See Appendix –Validated Advanced Sterilization Products (ASP) Cycles**

Gemini Wrap Weight	Advanced Sterilization Products STERRAD® Sterilization System Cycles
GEM0 - Super Light Weight	Not applicable
GEM1 - Light Weight	STERRAD® 50
GEM2 - Regular Weight	STERRAD® 100S
GEM3 - Medium Weight	STERRAD® 200S
GEM4 - Heavy Weight	STERRAD® NX™, (Standard Cycle, Advanced Cycle)
GEM5 - Super Heavy Weight	STERRAD® 100NX™, (Standard Cycle, Flex Cycle, EXPRESS Cycle, DUO Cycle)

**Table 6. STERIS V-PRO® Low Temperature Sterilization Systems.**  
**The wrap was validated to be effectively aerated during the pre-programmed cycles.**

Gemini Wrap Weight	STERIS V-PRO® Low Temperature Sterilization Cycles
GEM0 – Super Light Weight	Not applicable
GEM1 – Light Weight	STERIS V-PRO® 60 (Lumen, Non-Lumen and Flexible Cycles) STERIS V-PRO® 1 (Lumen Cycle) STERIS V-PRO® 1 Plus (Lumen and Non-Lumen Cycle) STERIS V-PRO® maX (Lumen, Non-Lumen, and Flexible Cycle)
GEM2 – Regular Weight	
GEM3 – Medium Weight	
GEM4 – Heavy Weight	
GEM5 – Super Heavy Weight	

**Table 7. STERIZONE® VP4 Sterilization System – See Appendix – Validated STERIZONE® VP4 Cycles**

Gemini Wrap Weight	STERIS V-PRO® Low Temperature Sterilization Cycles
GEM0 – Super Light Weight	Not applicable
GEM1 – Light Weight	STERIZONE® VP4 Sterilizer single preset sterilization cycle
GEM2 – Regular Weight	
GEM3 – Medium Weight	
GEM4 – Heavy Weight	
GEM5 – Super Heavy Weight	

The wrap is intended to allow sterilization of the enclosed medical device(s) and also maintain sterility of such content until used.

## WARNINGS

- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.
- Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient.
- Do not use wrap if damage or external particulates are identified prior to use.
- Do not use wrapped contents if the wrap is torn, wet, dirty or suspect of mishandling or abuse.
- Do not use Gemini Wrap in sterilization modes or conditions that fall outside the parameters listed in the Indications for Use section.

**Gemini Wrap is not intended for use in the following sterilization modes or cycles:**

- Dry heat sterilization
- Radiation (E.g., Gamma, E-Beam or other)
- Do not use wrap in sterilization cycles that will exceed 149°C.

## CAUTION

GEMO can be used only with the following sterilization methods : **Pre-vacuum steam**, and **Gravity steam**.

## PRECAUTIONS

The following are general precautions to be aware relating to Gemini Sterilization Wrap:

- Do not open case packaging with sharp knives to avoid creating accidental cuts in the wrap.
- Prior to being wrapped and sterilized, ensure that contents/devices have been deemed compatible for sterilization in the modes of sterilization identified in the above Indications for Use section. Consult the sterilization instructions for all devices intended for sterilization. Some medical devices, regardless of the sterilization method and sterilization wrap/container used, may require special consideration in packing configurations to ensure sterilization (refer to guidance EN ISO TS 16775 part 3.2.2.3 Assembly considerations).
- Do not use wrap in the presence of flammable materials. The Gemini wrap is non-conductive.
- If sterilization is performed by an outside contract facility, Medline Industries, Inc. recommends that the wrapped devices should be protected from contamination by an additional covering.

## INSTRUCTIONS FOR USE

Gemini Sterilization Wrap should be used in accordance with all preparation, wrapping, and sterilization chamber loading recommendations of the following standards:

- EN ISO 11607-1: Packaging for terminally sterilized medical devices; Requirements for materials, sterile barrier systems and packaging systems.
- EN 868-2 Packaging for terminally sterilized medical devices, Sterilization wrap, Requirements and test methods
- EN ISO TS 16775 Packaging for terminally sterilized medical devices – Guidance on the application of ISO 11607-1 and ISO 11607-2
- EN 285 Sterilisation – Steam sterilisers – Large sterilisers
- EN 13060 Small steam sterilisers
- EN 1422 Sterilisers for medical purposes, Ethylene oxide sterilisers, Requirements and test methods

## GENERAL STORAGE (PRE & POST STERILIZATION)

- Location should be clean, dust free and away from fluorescent or ultraviolet light.
- Use first in, first out (FIFO) stock rotation. The maximum sterility period has to be taken into account post sterilization (see page 8).
- Refer to guidance EN ISO TS 16775.

## PRIOR TO USE

- Examine wrap and discard if damage or extraneous matter is detected.
- Thoroughly clean and dry items to be wrapped/packaged.

# COMMON WRAPPING TECHNIQUES

- Place item(s) on wrap using typical aseptic wrapping techniques per EN ISO TS 16775 annex C. Other national standards or professional guidelines can be used for wrapping techniques. Recommendations for wrap contents are provided in Table 8.
- If using the simultaneous wrapping technique, ensure the first fold is pulled far enough to cover all package surfaces to ensure sterility maintenance.
- Secure the wrapped package with sterilization indicator tape or alternate closure method suitable for the sterilization method to be used.
- Closure must allow the sterilant to penetrate the wrapped package, avoid constriction of the package and maintain package integrity.

**Table 8: Wrap Model Recommendations<sup>1</sup>**

Gemini Wrap Weight	Intended Load	Maximum Recommended Wrapped Package Content <sup>2</sup>					
		Pre-Vacuum, Gravity	EO	ASP STERRAD® 50, 100S, 200S, NX™, and 100NX™ Cycles	STERIS V-PRO® 1, 1Plus, maX	STERIS V-PRO® 60	STERI-ZONE® VP4 Sterilizer
GEM0 – Super Light Weight	Not applicable	2.7 Kg	Not applicable				
GEM1 – Light Weight	Light weight package (for example: standard linen packs)	2.7 Kg	2.7 Kg	4.8 Kg	3 Kg	2.7 Kg	2.7 Kg
GEM2 – Regular Weight	Light to moderate weight package (for example: general use medical instruments)	4 Kg	4 Kg	4.8 Kg	4 Kg	4 Kg	4 Kg
GEM3 – Medium Weight	Moderate to heavy weight package (for example: general use medical instruments)	5.9 Kg	5.9 Kg	4.8 Kg	4.5 Kg	5.9 Kg	5.9 Kg
GEM4 – Heavy Weight	Heavy weight package (for example: general use medical instruments)	7.7 Kg	7.7 Kg	4.8 Kg	4.5 Kg	7.7 Kg	7.7 Kg
GEM5 – Super Heavy Weight	Very heavy weight package (for example: general use medical instruments)	11.3 Kg	11.3 Kg	4.8 Kg	4.5 Kg	11.3 Kg	11.3 Kg

<sup>1</sup> Individual results may differ due to factors such as variations in handling practices, wrapping techniques and folding methods. Results may also differ due to the use of irregularly shaped contents, which may put added stress on the wrap. Each healthcare facility should determine for itself which wrap model is most appropriate for each intended use.

<sup>2</sup> It is recommended to not exceed the maximum wrapped package content weights indicated for each wrap model. Furthermore is recommended to not exceed the number, weight and sizes of individual content types that were validated for the Gemini Sterilization Wrap.

## The following loads were used in the pre-vacuum steam and gravity Sterility Maintenance Validation Studies:

- Super Light Weight: 16 OR towels in two stack of eight
- Light Weight: 2 huck towels (43x73.5 cm), 3 fluid-resistant drapes (274.5x183 cm)
- Regular Weight: 16 huck towels (43x73.5 cm), 1 fluid-resistant table cover (228.5x152.5 cm). 2.3 Kg of metal mass.
- Medium Weight: 4 tray liners (51x63.5 cm) stacked 58.5x28x9 cm tray containing 5 Kg of metal mass.

- Heavy Weight: 4 tray liners (51x63.5 cm) stacked 58.5x28x9 cm tray containing 6.8 Kg of metal mass.
- Super Heavy Weight: 4 tray liners (51x63.5 in.) stacked 58.5x28x9 cm tray containing 9 Kg of metal mass.

#### **The following loads were used in EO Sterility Maintenance Validation Studies:**

- Super Light Weight: Not applicable
- Light Weight: 2 huck towels (43x73.5 cm), 3 fluid-resistant drapes (274.5x183 cm)
- Regular Weight: 16 huck towels (43x73.5 cm), 1 fluid-resistant table cover (228.5x152.5 cm). 2.3 Kg of metal mass.
- Medium Weight: 4 tray liners (51x63.5 cm) stacked 58.5x28x9 cm tray containing 5 Kg of metal mass.
- Heavy Weight: 4 tray liners (51x63.5 cm) stacked 58.5x28x9 cm tray containing 6.8 Kg of metal mass.
- Super Heavy Weight: 4 tray liners (51x63.5 in.) stacked 58.5x28x9 cm tray containing 9 Kg of metal mass.

#### **The following loads were used in the Advanced Sterilization Products (ASP) STERRAD® 50, 100S, 200S, NX™, and 100NX™ Sterility Maintenance Validation Studies:**

- Super Light Weight: Not applicable.
- Light Weight: 58.5x28x10 cm tray containing metal instruments.
- Regular Weight: 58.5x28x10 cm tray containing metal instruments.
- Medium Weight: 58.5x28x10 cm tray containing metal instruments.
- Heavy Weight: 58.5x28x10 cm tray containing metal instruments.
- Super Heavy Weight: 58.5x28x10 cm tray containing metal instruments.

#### **The following loads were used in the STERIS V-PRO® Sterility Maintenance Validation Studies:**

- Super Light Weight: Not applicable.
- Light Weight: 1.1 Kg metal mass, 6 forceps, STERIS V-PRO® Tray (43x25.5x9 cm) at 1.8 Kg.
- Regular Weight: 2.3 Kg metal mass, 6 forceps, STERIS V-PRO® Tray (43x25.5x9 cm) at 1.8 Kg.
- Medium Weight: 2.7 Kg. metal mass, 6 forceps, STERIS V-PRO® Tray (43x25.5x9 cm) at 1.8 Kg.
- Heavy Weight: 2.3 Kg metal mass, 6 forceps, STERIS V-PRO® Tray (53.5x25.5x9 cm) at 2.3 Kgs.
- Super Heavy Weight: 2.3 Kg metal mass, 6 forceps, STERIS V-PRO® Tray (53.5x25.5x9 cm) at 2.3 Kg.

#### **The following loads were used in STERIS V-PRO® 60 Sterility Maintenance Validation Studies:**

- Super Light Weight: Not applicable.
- Light Weight: V-PRO® Tray (25.5x53.5x9 cm), metal mass to make final total weight tested above, and 6 forceps.
- Regular Weight: V-PRO® Tray (25.5x53.5x9 cm), metal mass to make final total weight tested above, and 6 forceps.
- Medium Weight: V-PRO® Tray (25.5x53.5x9 cm), metal mass to make final total weight tested above, and 6 forceps.
- Heavy Weight: V-PRO® Tray (25.5x53.5x9 cm), metal mass to make final total weight tested above, and 6 forceps.
- Super Heavy Weight: V-PRO® Tray (25.5x53.5x9 cm), metal mass to make final total weight tested above, and 6 forceps.

#### **The follow loads were used in STERIZONE® VP4 Sterility Maintenance Validation Studies:**

- Super Light Weight: Not applicable.
- Light Weight: Symmetry Medical – MicroPack® Plastic Tray (38x25.5x4 cm), metal and non-metal instruments to make final total weight recommended above.
- Regular Weight: Symmetry Medical – MicroPack® Plastic Tray (38x25.5x4 cm), metal and non-metal instruments to make final total weight recommended above.
- Medium Weight: SteriPack Metal Tray (58.5x28x12.5 cm), metal and non-metal instruments to make final total weight tested recommended above.

- Heavy Weight: SteriPack Metal Tray (58.5x28x12.5 cm), metal and non-metal instruments to make final total weight tested recommended above.
- Super Heavy Weight: SteriPack Metal Tray (58.5x28x12.5 cm), metal and non-metal instruments to make final total weight tested recommended above.

## STERILIZATION PRACTICES

- Gemini Sterilization Wraps are intended for use with the common healthcare sterilization parameters listed in the Indications for Use. The sterilizer manufacturer should be consulted for appropriate sterilizer loading configurations.
- If a sterilizer malfunctions or a cycle is aborted before completion, packages should be re-wrapped prior to being placed into another sterilization cycle.
- Sets or items wrapped in Gemini Sterilization Wrap for sterilization should follow the load and stacking configurations identified in EN 285, EN 13060 or EN 1422 according to the sterilizer used.

NOTE: Many factors can affect drying time other than sterilization wrap, such as the pack configuration that is used, sterilizer loading, cycle variations, the performance of the sterilizer, temperature, altitude, and ambient environmental conditions. Sterilizers vary widely in design and performance characteristics. As recommended in the EN ISO TS 16775 guideline -part 3.3.2.9.3, the user should perform a qualification of the wrap and assess the contents after sterilization to ensure that sufficient drying has occurred in moist heat sterilization process.

## POST STERILIZATION COOLING/UNLOADING

- Leave wrapped packages on the sterilizer cart untouched until cool to avoid compromising package sterility.
- Visually inspect wrapped items as they are removed from the cart. Items that are torn, wet, or compressed should not be used.

## STERILITY MAINTENANCE

Under real time testing, Gemini Wrap is able to maintain package sterility for the following durations in controlled environmental conditions following sterilization. These durations do not prevent users from continuing to practice established healthcare facility protocols.

**Table 9: Gemini Wrap Sterility Maintenance**

Sterilization Modality	Maximum sterility period (until opened)
Pre-Vacuum & Gravity Steam	2 Years
Ethylene Oxide (ETO)	90 days
STERRAD® 50, 100S, 200S, NX™, and 100NX™ Cycles	1 Year
STERIS V-PRO® 1, 1Plus, maX, & 60	1 Year
STERIZONE® VP4	1 Year

# OPENING

- Inspect package for damage, moisture, or any sign of potential contamination prior to opening and again after opening but before use of contents. **Caution: Do not use contents if these conditions are present, as sterility could be compromised.** Reprocess the contents using an unprocessed wrap if any of these conditions are noted.
- Open packages aseptically in accordance with the healthcare facility's policy.

# DISPOSAL

- Do not re-use. Medline Industries, Inc. does not endorse the re-use (re-sterilization) of its sterilization wrap and does not warrant performance if product is re-used.
- Recycle, landfill or incinerate used wrap based upon state and local regulations. Only recycle non-soiled wraps.
- Note that the wrap is composed of polypropylene plastic which has a plastics recycling code of "5".
- Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

# APPENDIX:

**Note:** Refer to the User's Guide for complete instructions on load and cycle for each Sterilizer System below. The instructions provided below are not intended to replace the detailed Instructions for Use provide with each Sterilizer System.

## Validated Advanced Sterilization Products (ASP) STERRAD® 50, 100S, 200S, STERRAD® NX™, and STERRAD® 100NX™ Cycles

ASP STERRAD® System and Cycle	Intended Load
STERRAD® 50	<p>Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"><li>· An inside diameter of 1 mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens.</li><li>· An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens.</li><li>· An inside diameter of 3 mm or larger and a length of 400 mm or shorter of single-channel stainless steel lumens.</li><li>· An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel TEFLON®/Polyethylene lumens.</li></ul> <p>Refer to the STERRAD® 50 Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 10 lumens per load).</p>

ASP STERRAD® System and Cycle	Intended Load
STERRAD® 100S	<p>Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> <li>• An inside diameter of 1mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 3mm or larger and length of 400 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel TEFLON-Polyethylene lumens.</li> </ul> <p>Refer to the STERRAD® 100 Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 10 lumens per load).</p>
STERRAD® 200S	<p>Reusable metal and non-metal medical devices, including up to 12 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> <li>• An inside diameter of 1 mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 3 mm or larger and a length of 400 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel TEFLON®/Polyethylene lumens.</li> </ul> <p>Refer to the STERRAD® 200S Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 16.5 Kg per load).</p>
STERRAD® NX™ Standard Cycle	<p>Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> <li>• An inside diameter of 1mm or larger and a length of 150 mm or shorter of single-channel stainless steel lumens.</li> <li>• An inside diameter of 2 mm or larger and a length of 400 mm or shorter of single-channel stainless steel lumens.</li> </ul> <p>Refer to the STERRAD® NX™ Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 4.9 Kg per load).</p>

ASP STERRAD® System and Cycle	Intended Load
STERRAD® NX™ Advanced Cycle	<p>Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> <li>• An inside diameter of 1 mm or larger and length of 500 mm or shorter of single-channel stainless steel lumens.</li> </ul> <p>OR</p> <ul style="list-style-type: none"> <li>• One Single-Channel Flexible Endoscope with or without a silicone mat and no additional load. The flexible endoscope may contain:</li> <li>• A single-channel TEFLON®/Polyethylene lumen with an inside diameter of 1 mm or larger and a length of 850 mm or shorter.</li> </ul> <p>Refer to the STERRAD® NX™ Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 4.9 Kg per load).</p>
STERRAD® 100NX™ Standard Cycle	<p>Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> <li>• An inside diameter of 0.7 mm or larger and a length of 500 mm or shorter of single-channel stainless steel lumens. (A maximum of two flexible endoscopes, one per tray per sterilization cycle).</li> </ul> <p>Refer to the STERRAD® 100NX™ Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 9.7 Kg per load).</p>
STERRAD® 100NX™ Flex Cycle	<p>One or two single-channel Flexible Endoscope with or without a silicone mat and no additional load. The flexible endoscope may contain:</p> <ul style="list-style-type: none"> <li>• A single-channel TEFLON®/Polyethylene lumen with an inside diameter of 1 mm or larger and a length of 850 mm or shorter. (A maximum of two flexible endoscopes, one per tray per sterilization cycle).</li> </ul> <p>Refer to the STERRAD® 100NX™ Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 5.5 Kg per load).</p>
STERRAD® 100NX™ Express Cycle	<p>Non-lumened reusable metal and non-metal devices requiring surface sterilization, and sterilization of diffusion-restricted spaced such as the hinged portion of forceps and scissors, and rigid or semi-rigid endoscopes without lumens.</p> <p>Refer to the STERRAD® 100NX™ Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 4.9 Kg per load).</p>

ASP STERRAD® System and Cycle	Intended Load
STERRAD® 100NX™ DUO Cycle	<p>One or two single-channel Flexible Endoscope with accessory devices that are normally connected to it, with or without a silicone mat. The flexible endoscope may contain:</p> <ul style="list-style-type: none"> <li>• A single-channel TEFLON®/Polyethylene lumen with an inside diameter of 1 mm or larger and a length of 875 mm or shorter.</li> <li>• Accessory devices that are normally connected to a flexible endoscope during use.</li> <li>• Flexible endoscopes without lumens.</li> </ul> <p>Refer to the STERRAD® NX™ Sterilizer User's Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 4.9 Kg per load).</p>

## Validated STERIS V-PRO® Low Temperature Sterilization Systems Pre-programmed Cycle

STERIS® System and Cycle	Intended Load
STERIS® V-PRO® maX Lumen Cycle	<p>Reusable metal and non-metal medical devices including instruments with diffusion-restricted spaces (such as the hinged portion of forceps or scissors) and single, dual or triple channelled rigid/semi-rigid endoscopes, with the following configurations:</p> <ul style="list-style-type: none"> <li>• Single or dual channelled devices with stainless steel lumens with <ul style="list-style-type: none"> <li>◦ An inside diameter of 0.77 mm or larger and a length of 410 mm or shorter</li> </ul> </li> <li>• Triple Channelled devices with stainless steel lumens with <ul style="list-style-type: none"> <li>◦ An inside diameter of 1.2 mm or larger and a length of 257 mm or shorter</li> <li>◦ An inside diameter of 1.8 mm or larger and a length of 310 mm or shorter or</li> <li>◦ An inside diameter of 2.8 mm or larger and a length of 317 mm or shorter</li> </ul> </li> </ul>
STERIS® V-PRO® maX Flexible Cycle	<p>Single or dual channelled Flexible Surgical Endoscopes or Bronchoscopes with lumens that have:</p> <ul style="list-style-type: none"> <li>• An inside diameter of 1 mm or larger and a length of 990 mm or shorter.</li> </ul>
STERIS® V-PRO® maX Non-Lumen Cycle	<p>Reusable metal and non-metal non-lumened medical devices including non-lumened rigid, semi-rigid and flexible endoscopes and medical devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps or scissors.</p>

STERIS® System and Cycle	Intended Load
STERIS® V-PRO® 60 Lumen Cycle	<p>Reusable metal and non-metal medical devices including instruments with diffusion-restricted spaced (such as the hinged portion of forceps or scissors) and single, dual or trip channelled rigid/semi-rigid endoscopes, with the following configuration:</p> <ul style="list-style-type: none"> <li>Single or dual channelled devices with stainless steel lumens with <ul style="list-style-type: none"> <li>An inside diameter of 0.77 mm or larger and a length of 410 mm or shorter</li> </ul> </li> <li>Trip channelled devices with stainless steel lumens with <ul style="list-style-type: none"> <li>An inside diameter of 1.2 mm or larger and a length of 257 mm or shorter</li> <li>An inside diameter of 1.8 mm or larger and a length of 310 mm or shorter or</li> <li>An inside diameter of 2.8 mm or larger and a length of 317 mm or shorter</li> </ul> </li> </ul>
STERIS® V-PRO® 60 Flexible Cycle	<p>Single or dual channelled Flexible Surgical Endoscopes or Bronchoscopes with lumens that have:</p> <ul style="list-style-type: none"> <li>An inside diameter of 1 mm or larger and a length of 990 mm or shorter</li> </ul>
STERIS® V-PRO® 60 Non-Lumen Cycle	<p>Reusable metal and non-metal non-lumened medical devices including non-lumened rigid, semi-rigid, and flexible endoscope and medical devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps or scissors.</p>

### Validated STERIZONE® VP4 Sterilizer Single Preset Sterilization Cycle

STERIZONE® VP4 Sterilizer	Intended Load
STERIZONE® VP4 Sterilizer Single Preset Sterilization Cycle	<p>General medical instruments, representing the following geometries:</p> <ul style="list-style-type: none"> <li>Gliding mechanism, hinges, and screws, stopcock, instrument with rigid lumen (no dead end)</li> </ul> <p>Single channel flexible endoscopes with:</p> <ul style="list-style-type: none"> <li>Internal channel diameter of 1 mm or larger and length of 850 mm or shorter</li> </ul> <p>Single/Double Channel rigid endoscopes with:</p> <ul style="list-style-type: none"> <li>Internal channel diameter of 0.7 mm or larger and length of 500 mm o shorter</li> </ul> <p>Refer to STERIZONE® VP4 Sterilizer User's Guide for complete instructions on load(s), including chamber loading instructions.</p>

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Manufactured in a WRAP-certified facility.  
Medline supports labour and human rights.

Medical Device - Medizinprodukt - Dispositif médical - Medical device - Dispositivo Medico - Producto Sanitario -  
Dispositivo médico - Medicinteknisk Produkt - Lääkinnällinen laite - Medicinsk udstyr - Medisinsk utstyr -  
Zdravotnickým prostředkem - Dispozitiv medical - Orvostechnikai eszköz

Medline International France SAS,  
5 rue Charles Lindbergh,  
44110 Châteaubriant, France

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

Manufactured for:  Medline Industries, Inc., Three Lakes Drive, Northfield, IL 60093 USA. © 2021 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc.

RC21ENW

2021-03