en This product contains the phthalate DEHP (applicable for tubing only). Phthalates are widely used as plasticizers in medical products and have undergone extensive testing. However, when using the device, it should be taken into consideration that possible health effects when used for long durations in pregnant women, nursing mothers and children cannot be entirely discounted. Not Made With Natural Rubber Latex.

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

Medline Industries, Inc., Three Lakes Drive, Northfield, IL 60093 USA.

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia


Made in China. www.medline.com RB20JGP 2020-02

	1000		1500		3000
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	OR511 OR516 OR531
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	OR555 OR562 -
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	OR510 OR515 OR530
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	OR516K OR526K -
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	OR517K OR527K -
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	- OR529K -
	1000 cc/mL		1500 cc/mL		3000 cc/mL
				REF	- - OR539K


SETUP AND USE | SINGLE

A

B, C

D, E, F, G, H

en SEMI RIGID LINERS

READ INSTRUCTIONS COMPLETELY PRIOR TO USE.

Precautions:

- Avoid tilting, jarring, or sudden movements of the canister which may activate the hydrophobic filter and can cause premature shut-off.
- Allow for sufficient length of suction tubing (patient/liner) to prevent canister from tilting during use.
- For proper functioning of suction canisters, use appropriate vacuum tubing between canister and vacuum source that resist from kinking.
- Vacuum should only be introduced after correct device connections are made.
- Ensure vacuum suction remains on until all connecting tubing is disconnected and ports are securely capped.
- Always wear proper personal protective apparel when handling infectious waste.
- Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient.
- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

SETUP AND USE for Semi Rigid Liner : SINGLE

Port identification key

- ① Vacuum Port ② Patient Port ③ Ortho Port ④ Accessory Port

Directions for set up and use

- Remove liner from its packaging and snap lid and liner components together by using even pressure to fasten and make certain all tabs on lid are firmly snapped into place to ensure a secure fit.
- Insert proper sized semi rigid liner into the coinciding sized outer reusable canister. Connect tubing from the suction regulator (or other vacuum source) to the port underneath the canister.
- Press down on the liner lid in order to make sure that you have securely sealed the liner and outer canister.
- Attach red tubing from vacuum source to the Vacuum Port (1) on the semi rigid liner lid.
- Securely cap off Ortho Port (3) and Accessory Port (4)
- Attach patient tubing to the Patient Port (2) on the semi rigid liner lid.
- Confirm all unused ports are securely capped off prior to initiating suction.
- The system is ready for use.

Using the semi rigid liners

- When the suction is on, the liquid flows inside the semi rigid liner. The semi rigid liner can be used throughout the procedure or until the semi rigid liner is full. When the semi rigid liner is full, the shut-off valve cuts the suction and prevents the flow of liquid to the suction source.
- Each semi rigid liner is built with hydrophobic filter that blocks definitively when in contact with liquids. If canister is bumped or moved suddenly away from a vertical position, liquids might be in contact with the hydrophobic filter and create potential premature shut-off. So, if suction flow is stopped prematurely, change the semi rigid liner.

Using the semi rigid liners with solidifier

If adding solidifier, follow manufacturer's instructions. Solidifiers should be added before shutting off the vacuum source to prevent risk of fluid leaking from the accessory port when the solidifier is added.

As solidifier expands when absorbing fluids, measure accuracy is impacted.

Directions for disposal

Note: Vacuum should remain on until steps I-K are completed

- Disconnect the connective tubing from Patient Port (2).
- Ensure that all ports are capped off and sealed other than Vacuum Port (1).

K. Disconnect red vacuum tube from Vacuum Port (1) and cap that port securely.

L. Turn vacuum off.

M. Remove semi rigid liner from outer canister.

N. Comply with facility's protocols and local, state and federal requirements when disposing of semi rigid liners and contents.

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

⚠ WARNINGS

- Do not reuse or sterilize. If the same patient undergoes long-term treatment the semi rigid liner must be replaced at least every 24 hours.
- Do not subject the canister to active vacuum when not in use.
- Not intended for pleural/thoracic use.
- Inspect the device for the presence of cracks or damaged parts prior to use.
- Monitor vacuum gauge to assess for consistent vacuum pressure prior to and during use.
- Causes of premature shut-off include:
 - ◊ Bumping the canister or tilting at an angle can activate the shut-off mechanism inside the lid.
 - ◊ Suctioning agents that release gas, including hydrogen peroxide or povidone iodine may lead to the release of foam that can activate the the shut-off mechanism or even eject fluid from the canister.

Cleaning of reusable canister

Reusable canister can be wiped with disinfectant.

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Medline is a registered trademark of Medline Industries, Inc.