en This product contains the phthalate DEHP (applicable for tubing only). Phthalates are widely used as plasticizers in medical products and have undergone extensive testing. However, when using the device, it should be taken into consideration that possible health effects when used for long durations in pregnant women, nursing mothers and children cannot be entirely discounted. Not Made With Natural Rubber Latex.

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

Medline is a registered trademark of Medline Industries, Inc.
Manufactured for | fabriqué pour | fabricado para:
Medline Industries, Inc., Three Lakes Drive, Northfield, IL 60093 USA.
Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia
Made in China. www.medline.com RB20ZAT 2020-02

	1200 cc/mL 2000 cc/mL 3000 cc/mL	REF OR212 OR220 OR230
	1200 cc/mL 2000 cc/mL 3000 cc/mL	REF OR312 OR320 -
	1200 cc/mL 2000 cc/mL 3000 cc/mL	REF OR319 - -
	1200 cc/mL 2000 cc/mL 3000 cc/mL	REF OR395 - -
	1200 cc/mL 2000 cc/mL 3000 cc/mL	REF OR896 - -
	1200 cc/mL 2000 cc/mL 3000 cc/mL TANDEM	REF OR159 OR139 OR139 OR542

SETUP AND USE | SINGLE

SETUP AND USE | TANDEM

en RIGID CANISTER

READ INSTRUCTIONS COMPLETELY PRIOR TO USE.

Precautions:

- Avoid tilting, jarring, or sudden movements of the canister which may activate the hydrophobic filter and can cause premature shut-off.
- Allow for sufficient length of suction tubing (patient/liner) to prevent canister from tilting during use.
- For proper functioning of suction canisters, use appropriate vacuum tubing between canister and vacuum source that resist from kinking.
- Vacuum should only be introduced after correct device connections are made.
- Ensure vacuum suction remains on until all connecting tubing is disconnected and ports are securely capped.
- Always wear proper personal protective apparel when handling infectious waste.
- Reuse of single use devices creates a potential risk for the patient or user. It may cause contamination and/or impairment of function; which may lead to injury, illness and/or death of a patient.
- If a serious incident has occurred in relation to the device, it should be reported to Medline and the competent authority of the Member State in which the user and/or patient is established.

SETUP AND USE for Rigid Canister : SINGLE

Port identification key

- ① Vacuum Port
- ② Patient Port
- ③ Ortho Port
- ④ Accessory Port

Directions for set up and use

- A. Attached the blue lid to the top of the rigid canister by using even pressure to fasten and make certain all tabs on lid are firmly snapped into place to ensure a secure fit.
- B. Securely place rigid canister in the coinciding sized holder. Use the correlated sized canister ring for each sized canister. Connect tubing from the suction regulator (or other vacuum source) to the red manifold tubing underneath the holder.
- C. Attach red vacuum tubing on the holders to the vacuum port (1) on the rigid canister lid.
- D. Securely cap off Ortho port (3).
- E. Attach patient tubing to the Patient Port (2) on the rigid canister lid.
- F. Confirm all unused ports are securely capped off prior to initiating suction.
- G. The system is ready for use.

Using the rigid canister

- When the suction is on, the liquid flows into the rigid canister. The rigid canister can be used throughout the procedure or until the rigid canister is full. When the rigid canister is full, the shut-off valve cuts the suction and prevents the flow of liquid to the suction source.
- Each rigid canister is built with hydrophobic filter that blocks definitively when in contact with liquids. If canister is bumped or moved suddenly away from a vertical position, liquids might be in contact with the hydrophobic filter and create potential premature shut-off. So, if suction flow is stopped prematurely, change the rigid canister.

Using the rigid canister with solidifier

If adding solidifier, follow manufacturer's instructions. Solidifiers should be added before shutting off the vacuum source to prevent risk of fluid leaking from the accessory port when the solidifier is added. As solidifier expands when absorbing fluids, measure accuracy is impacted.

Directions for disposal

Note: Vacuum should remain on until steps H to J are completed

- H. Disconnect the connective tubing from Patient Port (2).
- I. Ensure that all ports are capped off and sealed other than Vacuum Port (1).
- J. Disconnect red vacuum tube from Vacuum Port (1) and cap that port securely.
- K. Turn vacuum off.
- L. Remove rigid canister from the holder.

Precautions should be taken when discarding this device and disposal of the device shall be made in accordance with applicable hospital or national regulations for biologically hazardous waste.

SETUP AND USE FOR RIGID CANISTERS : TANDEM

Port identification

- ① Vacuum Port
- ② Patient Port
- ③ Ortho Port
- ④ Accessory Port

Directions for set up and use

- In a tandem application canister orientation should follow a clockwise pattern.
- A. Place rigid canisters in holders on the canister carrier.
- B. Ensure that the red manifold tubing is attached to each of the canister holders.
- C. Connect tubing from the suction regulator (or other vacuum source) to the red manifold tubing underneath the holder.
- D. Connect tandem tubing from Ortho Port (3) of the first canister to Patient Port (2) of second canister (removing elbow connector) that the overflow will enter. In order to ensure a secure seal, press down firmly.
- E. Repeat the tandem tube instruction above (Step D) until all the canisters are connected.
- F. Attach connecting patient suction tubing to the Patient Port (2) of first canister and attach adapted suction device (eg. Yankauer) to other end of the patient suction tubing.
- G. Confirm all unused ports are securely capped off prior to initiating suction.
- H. The system is ready for use.

Using the rigid canister with solidifier

Follow solidifier manufacturer's instructions for use. It is not recommended to use solidifier prior to suction during tandem configuration, as fluid will solidify, volume expands and suction is altered in the tandem tube during procedures. Solidifiers should be added before shutting off the vacuum source to prevent risk of fluid leaking from the accessory port when the solidifier is added.

Directions for disposal

Note: Vacuum should remain on until steps I to L are completed

- I. Disconnect the patient tubing from the Patient Port (2) on the first rigid canister (removing elbow connector) and securely cap the port.
- J. Carefully remove all tandem tubes between canisters and cap ports.
- K. Disconnect the red vacuum tubes coming from the holder from the Vacuum Port (1) and cap the port of the lid. Repeat this step until all red vacuum tubes are disconnected.
- L. Confirm that all ports are securely capped off.
- M. Turn off the vacuum.
- N. Remove the canisters from the holder carefully.
- O. Comply with facility's protocols and local, state and federal requirements when disposing of rigid canisters and contents.

Note: canister contents are considered hazardous and proper caution should be taken when handling. For rigid canisters tandem setups, it is recommended that all used and unused disposable canisters be discarded and subsequently replaced following each procedure as they are all contaminated by aerosols.

! WARNINGS

- Do not reuse or sterilize. If the same patient undergoes long-term treatment the rigid canister must be replaced at least every 24 hours.
- Do not subject the canister to active vacuum when not in use.
- Not intended for pleural/thoracic use.
- Inspect the device for the presence of cracks or damaged parts prior to use.
- Monitor vacuum gauge to assess for consistent vacuum pressure prior to and during use.
- Causes of premature shut-off include:
 - ◊ Bumping the canister or tilting at an angle can activate the hydrophobic filter inside the lid.
 - ◊ Suctioning agents that release gas, including hydrogen peroxide or povidone iodine may lead to the release of foam that can activate the hydrophobic filter or even eject fluid from the canister.

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