
## MED-GEL

en	INSTRUCTIONS FOR USE
de	GEBRAUCHSANWEISUNG
fr	NOTICE D'UTILISATION
nl	GEBRUIKSAANWIJZING
it	ISTRUZIONI PER L'USO
es	INSTRUCCIONES DE USO
pt	INSTRUÇÕES PARA USO
sv	INSTRUKTIONER FÖR ANVÄNDNING
fi	KÄYTTÖOHJEET
da	INSTRUKTION
no	BRUKSANVISNING

15 g

2 Kg

 **Medline International France SAS,**  
5 rue Charles Lindbergh, 44110 Châteaubriant, France.

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia  
RJ19HBN 2019-10

**Description:**

MED-GEL is a granular substance (Sodium polyacrylate, cross linked) which jellifies aspirated fluids in the Operating Theatre and Intensive Care Unit to facilitate safe disposal of waste.

**Composition/Ingredient:**

Sodium Polyacrylate (95%), Water (5%).

**Recommended Dosage:**

RE2015: MED-GEL is packed in an individual 15gr water-soluble bag. 1 bag = 1 dose = 15gr MED-GEL

RE2050: MED-GEL is packed in a 2kg container which holds 2 x1kg bag and a measuring cup. Use the measuring cup to pour MED-GEL to jellify the fluids. 1 full measuring cup = 1 dose = 15gr MED-GEL

Volume of fluids to be jellified	Quantity of MED-GEL required	Equivalence RE2015	Equivalence RE2050
1.0 L	15 gr of MED-GEL	1 dose of 15 gr	1 cup of 15 gr
1.5 L	30 gr of MED-GEL	2 doses of 15 gr	2 cups of 15 gr
3.0 L	45 gr of MED-GEL	3 doses of 15 gr	3 cups of 15 gr

The fluids are jellified after approximately 10 - 15 minutes.

Dosage is dependent on the composition of the fluids. A high percentage of electrolytes in the solution typically requires a higher dosage. Gallbladder fluids can partially reduce the effectiveness (fat).

**Recommended Application: with disposable suction liner and canister system**

MED-GEL is added in the appropriate quantity to the fluids in the disposable suction liner or canister. It is recommended to add MED-GEL (via the accessory port) into the suction liner after suctioning and prior switching off the vacuum source.

- In single configuration: MED-GEL may be added into the liner prior to suctioning. If premature shutoff occurs it will not be reversible and will require the immediate change of the liner.
- In tandem configuration: Do not put MED-GEL into the liner prior to suctioning since this may reduce suction effectiveness. Add MED-GEL after suctioning and prior switching off the vacuum source.

**First aid measures:**

- Skin Exposure: In case of contact, immediately wash skin with soap and copious amounts of water. If irritation persists, call a physician.
- Eye exposure: In case of contact with eyes, flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating the eyelids with fingers. If irritation persists, call a physician.
- Inhalation Exposure: If inhaled, away from the area of exposure. If necessary, seek medical attention.
- Oral Exposure: If swallowed, irrigate the mouth with water and seek medical attention.

**Exposure controls/Personal protection:**

- Respiratory protection: Use respiratory protection when dust is formed. Ensure sufficient workplace ventilation.
- Eye protection: Protective goggles
- Hand protection: Not applicable

**Fire-fighting measures:**

- Extinguishing media: Suitable: water spray, dry chemical, carbon dioxide or appropriate foam.
- Firefighting protective equipment: wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.

**Toxicity, Dangers:**

- Acute oral toxicity: LD50, dose > 5000 mg/kg, limit test.
- Skin irritation: No irritation.
- Eye irritation: Very slight eye irritation. Particle effect.
- Sensitization: No sensitization.
- Acute dermal toxicity: LD50, dose > 2000 mg/kg, limit test.

**Disposable Considerations:**

- To be stored under dry conditions.
- Do not dispose of jellified fluids via the sewage system. Comply with facility's protocols and local, state and federal requirements for disposal (COSHH regulations for UK only).
- In regards to the present knowledge the product is a non-dangerous substance in the definition of the European Regulation N°1272/2008 – Classification, Labeling and Packaging of substances and mixtures.

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