




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



Technical sheet

Back table cover

Reference(s)

1900CEB, 2100CEB, 8186CEB, 8368CEA, 8371CEB, 8378CEB, 8380CEB, DYJPETC001, DYJPETC060, DYJPETC065, DYNJPE2309, DYNJPE2312, DYNJPE2314, DYNJPE2317, DYNJPE2320

Description

A SMS wrap contains :

Reference	Width « A » (cm)	Height « B » (cm)	PE Thickness (µm)	Reinforcement	Qty per carton
1900CEB	140	190	75	SPP reinforcement 74x190cm	22
2100CEB	140	120	75	SPP reinforcement 66x120cm	26
8186CEB	160	240	90	SPP reinforcement 74x240cm	14
8368CEA	140	120	50	SPP reinforcement 64x120 cm	30
8371CEB	140	190	50	SPP reinforcement 74x190cm	22
8378CEB	125	230	50	SPP reinforcement 74x230cm	18
8380CEB	203	240	75	SPP reinforcement 100x240cm	14
DYJPETC001	110	112	50	SPP reinforcement 64x112 cm	20
DYJPETC060	152	190	50	SPP reinforcement 74x190 cm	22
DYJPETC065	152	152	50	SPP reinforcement 74x152 cm	22
DYNJPE2309	110	190	50	SPP reinforcement 64x190 cm	22
DYNJPE2312	196	229	60	Wetlaid reinforcement 102x229 cm	22
DYNJPE2314	152	229	90	Wetlaid reinforcement 76x229 cm	12
DYNJPE2317	201	279	90	Polyrayon reinforcement 102x279 cm	12
DYNJPE2320	152	230	50	SPP reinforcement 74x230 cm	22

Main materials

Blue PE (polyethylene) with reinforced area:
 - Polyrayon absorbent reinforcement
 - Repellent SMS (polypropylene)
 (polyester + viscose)
 - Wetlaid (cellulose)
 - Absorbent SPP (polypropylene)

Packaging

Peel pouch or Header bag

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
 In compliance with the standard EN13795
 MDD Class: I sterile



STERILE EO



© 2016 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_TableCover01_EN02

01.01.INT.12