




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



ECLIPSE®

Lithotomy pack

Reference(s)

29014CEA

Description

A sterile wrap 86x86 cm contains:

- 1 back table cover, 140x190 cm, 50 µm containing:

- 1 spunlace lithotomy drape 160/89x300 cm, 48x49 cm reinforcement, 6x17 cm fenestration with adhesive, 2 leggings, fluid control pouch, 2 cord holding tabs
- 1 absorbent towel, 30x20 cm unfolded

Main materials

Eclipse® spunlace (polyester + woodpulp)

Polyrayon absorbent reinforcement (polyester + viscose)

Pack factor

7 units per carton

Packaging

Vented sterilization pouch

1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS

EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC

In compliance with the standard EN13795

MDD Class: I sterile



Latex free



0086



STERILE|EO

© 2013 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-29014CEA-EN03

01.01.INT.12