




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



Technical sheet

Sleeves with cuffs

Reference(s)

598CE

Description

Breathable and impervious sleeves, with cuff, 14x62cm

Main materials

Prevention Plus 78gsm microporous and breathable

Tested against the following chemotherapy drugs per ASTM F739-99a:

- Dacarbazine
- Carmustine (BiCNU)
- Cyclophosphamide (Cytosan)
- Doxorubicin Hydrochloride
- 5- Fluorouracil
- Cisplatin
- Etoposide (Toposar)
- Paclitaxel (Taxol)
- Thio- Tapa
- Mitomycin C
- Vincristine Sulfate

Pack factor

80 units per carton

Packaging

Paper and film peel pouch – 2 pieces per pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile



Latex free



0086



STERILE EO

© 2016 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_GenSurgDrape_598CE_EN02

01.01.INT.12