

REF 9716CE


x24

Eclipse®

(01) 50 88438916388 5(17) YYMM28(10) YYMRPxxx

- en Surgical Gown
- fr Casaque de Chirurgie
- de OP-Mantel
- it Camice Chirurgico
- es Bata Quirúrgica
- pt Bata Cirúrgica
- nl Operatiejas
- sv Operationsrock
- da Operationskittel
- fi Leikkaustakki
- el Χειρουργική Μπλούζα
- no Operasjonsfrakk
- hu Operációs köpeny
- ro Halat chirurgical
- tr Cerrahi Önlük
- pl Fartuch chirurgiczny
- cs Operační plášť
- ru Хирургический халат

XLL
150 cm



LOT



0086

STERILE EO



LAB161305
Made in China

Medline International France SAS
5 rue Charles Lindbergh, 44110 Châteaubriant, France
Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

Label information (8"x6"): YYMM28: expiration date
5 years guarantee, expiration date displays
following date: 4 years, 11 months, end of month

YYMRPxxx: lot number
YY for Year of Manufacturing
M for Month of Production
R for Product Specification revision
P for Plant of Manufacturing
xxx for order incrementation

color: coated label, purple PMS2597


barcode created according to GS1-128 EAN



REF 9716CE

Eclipse®

- en Surgical Gown
- fr Casaque de Chirurgie
- de OP-Mantel
- it Camice Chirurgico
- es Bata Quirúrgica
- pt Bata Cirúrgica
- nl Operatiejas
- sv Operationsrock
- da Operationskittel
- fi Leikkaustakki
- el Χειρουργική Μπλούζα
- no Operasjonsfrakk
- hu Operációs köpeny
- ro Halat chirurgical
- tr Cerrahi Önlük
- pl Fartuch chirurgiczny
- cs Operační plášť
- ru Хирургический халат



STERILE/EO



LOT



(01) 00 88438916388 0(17) YYMM28(10) YYMRPxxx

LAB161306



Medline International France SAS

5 rue Charles Lindbergh, 44110 Châteaubriant, France

Distributed in Australia by: Medline International Two Australia Pty. Ltd., Sydney, Australia

Label information (14x21cm):

YYMM28: expiration date
5 years guarantee, expiration date displays
following date: 4 years, 11 months, end of month

YYMRPxxx: lot number
YY for Year of Manufacturing
M for Month of Production
R for Product Specification revision
P for Plant of Manufacturing
xxx for order incrementation

color: uncoated label, purple PMS2597

IFU to be used: IFU_TES02,
printed at the back of the insert label

Traceability stickers: applied on
the back of the ste pouch

barcode created according to GS1-128 EAN

Barcode should be printed in black

REF 9716CE

www.medline.com

LOT



REF 9716CE

www.medline.com

LOT



REF 9716CE

www.medline.com

LOT



(01)00884389163880

REF 9716CE

www.medline.com

LOT



(01)00884389163880