





European-Made Facemasks





Our surgical and procedural facemasks range is designed to provide high levels of safety and comfort. Tailored with healthcare professional and everyday use in mind, these masks meet rigorous European standards.

Procedural facemasks




| | Item number | Bacterial filtration efficiency (%) | Differential pressure (Pa/cm ²) | Splash resistance (kPa) | Packaging |
|---|---|-------------------------------------|---|-------------------------|----------------|
|  | BFE-230 | 99.8 | 49.6 | | 50/bx, 2000/cs |
| | Ideal for users with sensitive skin Spunbond polypropylene interior and exterior, white, with ear loops | | | | Type II |
|  | BFE-231 | 99.8 | 49.6 | | 50/bx, 2000/cs |
| | Spunbond polypropylene interior, spunbond polypropylene exterior, blue, with ear loops | | | | Type II |
|  | BFE-245 | 99.8 | 49.6 | 16 | 50/bx, 2000/cs |
| | Spunbond polypropylene interior and exterior, purple, with ear loops | | | | Type IIR |
|  | BFE-240 | 99.8 | 49.6 | 16 | 50/bx, 2000/cs |
| | Ideal for users with sensitive skin Spunbond polypropylene interior, spunbond polypropylene exterior, white, ear loops | | | | Type IIR |

Surgical facemasks

| | Item number | Bacterial filtration efficiency (%) | Differential pressure (Pa/cm ²) | Splash resistance (kPa) | Packaging |
|---|---|-------------------------------------|---|-------------------------|----------------|
|  | BFE-345 | 99.8 | 49.6 | 16 | 50/bx, 2000/cs |
| | Spunbond polypropylene interior, spunbond polypropylene exterior, purple, with ties | | | | Type IIR |
|  | BFE-346 | 99.8 | 49.6 | 16 | 50/bx, 2000/cs |
| | Spunbond polypropylene interior, spunbond polypropylene exterior, pink, with ties | | | | Type IIR |

Respirators

Respirators have been developed to provide high-level protection for healthcare professionals and patients by helping to reduce exposure to airborne pathogens and providing comfort during long shifts. These masks are available in two or three folds and with different protection levels (FFP2 or FFP3) to fit various needs.


| | Item number | Total inward leakage (%) | Sodium chloride test 95 l/min., % max. | Paraffin oil test 95 l/min., % max. | Packaging |
|--|--|--------------------------|--|-------------------------------------|---------------|
|  | 1220 | 0.85 | 0.47 | 1.61 | 50/bx, 200/cs |
| | Flat-folded, 100% polypropylene interior and exterior, dual polypropylene/elastane headbands, white, single-packed | | | | FFP2 |
|  | 2120L | 1.95 | 0.4 | 0.57 | 25/bx, 200/cs |
| | Tri-folded, 100% polypropylene interior and exterior, dual polypropylene/elastane headbands, white, single-packed | | | | FFP2 |
|  | 2130 | 1.95 | 0.4 | 0.57 | 25/bx, 200/cs |
| | 100% polypropylene interior and exterior, dual polypropylene/elastane headbands, white, tri-folded and single-packed | | | | FFP3 |



Experience excellent protection and better lead time, with our affordable **European-made facemasks.**



Raw materials sourced from **European countries** enables us to be reliable but flexible with delivery.



For more information, please contact your Medline representative or visit our website www.medline.eu/uk

Procedural Facemasks and Surgical Facemasks are class I non-sterile medical devices.

These products are intended for use by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2025 Medline Industries, LP. Medline is a registered trademark of Medline Industries, LP. Three Lakes Drive, Northfield, IL 60093, USA. ML1575_EN 04/2025

INSPIRE
Imp Helène Boucher, ZA Les Béluzes
42720 Pouilly /s Charlieu France

