



CENTURION
MEDICAL PRODUCTS


Compass UniversalHG for IAP Monitoring

Monitor intra-abdominal pressure quickly and accurately

**Precise
pressure
made simple**



IAP monitoring



The reliable measurement of bladder pressure is crucial to early detection, treatment and prevention of intra-abdominal hypertension (IAH) and abdominal compartment syndrome (ACS). Compass enables clinicians to easily and efficiently measure intra-abdominal pressure (IAP) with confidence.

Benefits and features of Compass

- **Intuitive/easy to use:** Replaces traditional pressure transducers; LCD screen displays easy-to-read digital pressure readings
- **Minimises risk of infection:** Attaches directly to the sampling port to keep the system closed and promote sterility
- **Accurate, reproducible measurements per ANSI standards:** Assists with delivering timely and precise IAH or ACS diagnoses for pressures ranging from -199 mmHg to +999 mmHg
- **Fast reading:** In as little as 1 – 3 seconds, saving time and promoting faster treatment

Ordering information

Item number	Description	Packaging
CUHG1-INT	Compass UniversalHG	5 ea/cs

For more information on this product, please contact your Medline account manager or visit our website: www.medline.eu/uk



Medline Industries Ltd
3rd Floor
Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom
Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk | uk-customerservice@medline.com

FOLLOW US  

Compass UniversalHG is a class IIa medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling.



All rights reserved. We reserve the right to correct errors that may have occurred.

© 2020 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. ML1015-EN 05/2020