

Pulse Lavage

For wound cleaning and debridement




Medline's pre-assembled Pulse Lavage System is a ready-to-use solution, designed for use in several applications, including orthopaedic, trauma surgery preparation and wound hydrotherapy.




Cone tip
Coaxial tip



All components are contained within a single, sterile blister pack, prepared for instant use.



Integrated suction and irrigation with suction tubing, designed to be compatible with all suction systems.



Integrated recyclable external alkaline battery that is easy to discard.

Ordering information

Item number	Description	Packaging
201-04-00	Pulse Lavage System	6 ea/cs

For more information, please contact your Medline account manager, or visit our website: www.medline.eu/uk

Pulse Lavage System is a class IIa sterile medical device intended for use by healthcare professionals. Before use, consult the instructions and precautions on the corresponding labelling. All rights reserved. We reserve the right to correct errors that may have occurred.



WELLKANG LTD.
Enterprise Hub, NW Business
Complex, 1 Beraghmore Rd,
Derry, BT48 8SE, Northern Ireland

