

Economy Sahara OR Table Sheet

OR set-up made easier




Medline's QuickSuite® OR Table Kits contain an assortment of disposable components for use before and after procedures to set-up and rapidly clear the OR table, respectively.

Our Economy Sahara Sheet is one of the components that help you conquer some of your largest OR challenges by:

- Boosting efficiency
- Promoting safety
- Providing a cost-effective solution

Key features

- Core polymer technology that wicks moisture away, maintaining skin integrity
- One size for OR table use: 102 x 229 cm
- Quilted, smooth surface
- Disposable design creates low risk for harbouring bacteria
- Glued edges help avoid linting
- Impervious back sheet



Ordering information

Item number	Description	Quantity
DYND4090E	Economy Sahara Sheet	20/case

For more information on this product, please contact your Medline account manager or visit our website: www.medline.eu/uk

Economy Sahara sheet is a class I NS medical device intended to be used by healthcare professionals. Before use, consult instructions and precautions on the corresponding labelling. We reserve the right to correct errors that may occur within this brochure.

© 2020 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. ML566-EN 05/2020.

Medline International France SAS
5, rue Charles Lindbergh
44110 Châteaubriant, France

