

Urinary Catheterization Kit

Range UROKIT00x



UROKIT001

- 1 drape 75 x 90cm
- 1 drape 75 x 90 cm with a 9 x 12 cm fenestrated
- 1 plastic tray
- 5 non-woven swabs, 7.5 x 7.5 cm
- 1 pair of Kocher forceps

Components basic supplements in the kit:

- 1 lubricating Jelly 5 mg
- 1 pre-filled syringe with 20 ml of water



UROKIT001-1

- 1 drape 75 x 90cm 40gsm
- 1 drape 75 x 90 cm fenestrated 9 x 12 cm adhesive and pre-cut 40gsm
- 1 plastic tray
- 5 non-woven swabs, 7.5 x 7.5 cm
- 1 pair of Kocher forceps
- 1 lubricating Jelly 5 mg

Components basic supplements in the kit:

- 1 pre-filled syringe with 20 ml of water



UROKIT002

- 1 rigid blister pack
- 1 drape 75 x 90cm 40gsm
- 1 drape 75 x 90 cm fenestrated 9 x 12 cm adhesive and pre-cut 40gsm
- 5 non-woven swabs, 7.5 x 7.5 cm
- 1 pair of Kocher forceps
- 1 lubricating Jelly 5 mg

Components basic supplements in the kit:

- 1 pre-filled syringe with 20 ml of water



UROKIT003 – Urinary Catheterization Kit with Preparation and Cleansing Component

Components for cleansing

In ablistar with 1 compartment:

- 1 non woven wash glove
- 2 exam gloves (size 7/8)
- 10 non-woven swabs 7.5 x 7.5 cm
- 2 absorbent hand towel 40 x 60 cm

Components for catheterization

- 1 drape 75 x 90 cm
- 1 drape 75 x 90 cm fenestrated 9 x 12 cm adhesive and pre-cut 40gsm
- 1 tray with 2 compartments
- 1 pair of Kocher forceps
- 10 non-woven swabs 7.5 x 7.5 cm
- 1 absorbent hand towel 40 x 60 cm
- 2 adhesive strips 14.8 x 1.25 cm
- 1 lubricating Jelly 5 mg

Components basic supplements in the kit:

- 1 pre-filled syringe with 20 ml of water



UROKIT004

- 1 rigid blister
- 1 drape 75 x 90cm 40gsm
- 1 drape 75 x 90 cm fenestrated 9 x 12 cm adhesive and pre-cut 40gsm
- 10 non-woven swab, 7.5 x 7.5 cm
- 1 pair of Kocher forceps
- 1 lubricating Jelly 5 mg

Components additional in the kit:

- 1 syringe pre-filled with 20 ml water



For more information, please contact your Medline representative or visit our website www.medline.eu/uk

Before use, please consult instructions and precautions on the corresponding labelling.

All rights reserved. We reserve the right to correct errors that may have occurred.

© 2026 Medline. ML1781_EN 08/2026