



Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



Technical sheet

Bouffant cap

Reference(s)

Reference	Color	Size	Main materials	Basis weight
CRIE1000	Green	Medium ± 53 cm	Spunbond polypropylene	12 gsm
CRIE1002	White			
CRIE1003	White	Large ± 60 cm	SMS (polypropylene)	10 gsm
CRIE1004	Blue			
CRIE1005	Green			
NONE28233F	Butterfly print			
NONE28234P	Multi-color print		Spunbond polypropylene	20 gsm

Description

Bouffant cap

Pack factor

500 units per carton

Packaging

5 dispenser boxes contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: 1 Non-Sterile



Latex free



© 2016 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_Headwear01_EN01 01.01.INT.12