



Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



CRIE2010

NONE28858 / NONE28859

CRIE2002

Technical sheet

Shoe covers

Reference(s)

Reference	Description	Color	Dimensions	Main materials	Basis weight	Packaging
CRIE2002	Shoe covers with textured sole	Blue	L 41 cm	Spunbond polypropylene	35g/m ²	100/bx 300/cs
CRIE2010	Shoe covers	Blue	One Size Fits Most 40 cm	Polyethylene	27g/m ²	10/roll 100/bag 1000/cs
NONE28858	Shoe covers with textured sole	Blue	L 40 cm	SMS (polypropylene)	35g/m ²	100/bx 300/cs
NONE28859			XL 43 cm			100/bx 200/cs

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical device regulation 2017/745
MDD Class: 1 Non-Sterile

