



NON-STERILE ISOLATION GOWNS

TECHNICAL DATA



Product No.	Size	Weight	Length	Packaging	Quantity
CRIE4000	Regular/Large	20 grams SPP	114 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE4001	X-Large	20 grams SPP	127 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE4000B	Regular/Large	20 grams SPP	114 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE4000G	Regular/Large	20 grams SPP	114 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE4001CF	X-Large	20 grams SPP	127 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE4010CF	Regular/Large	40 grams PE coated PP	114 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE4011CF	X-Large	40 grams PE coated PP	127 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
CRIE5000	Regular/Large	25 grams PE	114 cm <i>(minimum)</i>	Individual Poly Bag, Inner Box, Corrugate Case	15/bx, 75/cs
CRIE5001	X-Large	25 grams PE	134 cm <i>(minimum)</i>	Individual Poly Bag, Inner Box, Corrugate Case	15/bx, 75/cs
NONE27SMS3	Regular/Large	25 grams SMS	114 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	100/cs
NONE27SMS3XL	X-Large	25 grams SMS	127 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	100/cs
NONE27116	Regular/Large	48 grams PE coated PP	114 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
NONE27116XL	X-Large	48 grams PE coated PP	127 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs
NONE27279OHL	X-Large	48 grams PE coated PP	127 cm <i>(minimum)</i>	Inner Poly Bag Corrugate Case	50/cs

Description:	Gown protects healthcare worker and patient
Material:	Gowns made from: 100% Spunbond Polypropylene Nonwoven Material, 100% SMS Polypropylene Nonwoven Material, Polyethylene (PE) Material or Polyethylene Coated Polypropylene Material
Storage Conditions:	At room temperature
Packaging:	Ten gowns per poly bag (unsealed), 5 or 10 bags per case or 15 gowns per inner box, 5 boxes per case
Sterilization:	Non-Sterile
Manufacturer:	Medline Industries, Inc Conforms to ISO 13485
Regulation:	Conforms to Medical Devices Directive 93/42/EEC - PPE 89/686/EEC Class: 1



Latex free

Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.

Medline Industries, Inc
The Company of Choice
www.medline.com

 Medline
One Medline Place
Mundelein, Illinois 60060 USA
Telephone +1 800 MEDLINE
International@medline.com

LITIsolationGowns-EN01