




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



Technical sheet

Arthroscopy drape

Reference(s)

DYJPEARDTB

Description

A sterile wrap, 112x112 cm, contains:

- 1 tri-laminate arthroscopy drape, 231/290x315 cm, 6 cm diameter elastic fenestration, fluid collection pouch with exit ports, cord holding tab, armboard covers

Main materials

Absorbent and impervious tri-laminated material (polyethylene + polypropylene)

Pack factor

6 units per carton

Packaging

Vented sterilization pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile



Latex free



0123



STERILE EO



© 2013 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-DYJPEARDTB-EN01

01.01.INT.12