




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ ADVANCED

Double split sheet

Reference(s)

DYJPEDSDSM

Description

A SMS wrap, 91x91 cm contains :

- 1 SMMMS double split drape 240x286 cm, 2 adhesives split 6x90 cm, 2 absorbent reinforcements 35x68 cm and 1 absorbent reinforcement 70x89 cm, 2 hook & loop tube holders

Main materials

Repellent SMMMS (polypropylene)

Polyrayon absorbent reinforcement (polyester + viscose)

Packaging

Vented sterilization pouch

8 units per case

Manufacturer

Medline Industries Inc Mundelein USA

EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: I sterile



Latex free



0086



STERILE EO

© 2013 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-DYJPEDSDSM-
EN01

01.01.INT.12