



**MASTER
CARTON
MUST BE
OPENED**

Each Item:

REF DYND70802UK

Contents: 20 eaches per Master Carton

PO:



LOT

CE **EC REP** Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

Manufactured in China for: **Medline Industries, Inc.**,
Three Lakes Drive, Northfield, IL 60093 USA RJ18WZA

rev 2/11



**MASTER
CARTON
MUST BE
OPENED**

Each Item:

REF DYND70802UK

Contents: 20 eaches per Master Carton

PO:



LOT

CE **EC REP** Medline International France SAS, 5 rue Charles Lindbergh, 44110 Châteaubriant, France

Manufactured in China for: **Medline Industries, Inc.**,
Three Lakes Drive, Northfield, IL 60093 USA

rev 2/11

File Name: PKG_DYND70800UK_MSC_WZA.ai **PROD MGR** _____

Designer: Caroline

Date Started: 24/10/18

Date Finished:



QA _____

APPROVED APPD W/CHANGES REPROOF

PMS COLORS USED

CMYK

Revision: 0

Black

199



Note for the printing:

UDI area, generate and apply at manufacture

LOT number format: If possible, please follow below rule.

XXXYYWWN: Vendor code (XXX) - Year (YY) - W (Week) - WN (for week number)

If Gross Weight, Net Weight and PO are not on the printing files, those are to be added somewhere on the case (printed/stamped/sticked...) Please don't forget to add them at the time of physical proof/manufacture



Clipper base

Clipper base



REF DYND70802UK

ELECTRIC SURGICAL
CLIPPER, BASE



Medical Device

Manufactured in China for:
Medline Industries, Inc.,
Three Lakes Drive, Northfield, IL 60093 USA

EC REP

Medline International France SAS,
5 rue Charles Lindbergh,
44110 Châteaubriant, France

RJ18WZA

Date of Mfg (11): YYYY-MM-DD

LOT Lot (10): XXXYYWWN

P.O.:



(01) 10884389962954



(11)YYMMDD(10)XXXYYWWN

File Name: PKG_DYND70802UK_ECH_WZA.ai PRD MGR _____

Designer: Caroline

Date Started: 24/10/18

Date Finished:



QA _____

APPROVED APPD W/CHANGES REPROOF

PMS COLORS USED

CMYK

Revision: 0

Black



Note for the printing:

UDI area, generate and apply at manufacture

LOT number format: If possible, please follow below rule.

XXXYYWWN: Vendor code (XXX) - Year (YY) - W (Week) - WN (for week number)

4 x 2 cm

Charge Medline Specified Surgical Clipper only

MEDLINE Surgical Clipper Charger

Model - DYND70802UK

Input - DC 8V 200mA

Output - DC 3.3V 100mA

Conforms to: IEC 60601-1

Made in China

RI17WZA

LOT



ETL CLASSIFIED



Intertek



PMS COLORS USED



Revision: 0

Black

