



Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



Technical sheet

Prevention™ Plus Surgical Gown

Reference(s)	Reference	Size
	DYNJPE2301P	L, 114 cm
	DYNJPE2302P	XLL, 140 cm
Description	A SMS wrap 61x61 cm contains: - 1 hand towel 69x37 cm, unfolded - 1 Prevention™ Plus surgical gown, hook and loop adjustable neck	
Main materials	Prevention™ Plus trilayer fabric (spunbond polypropylene + polyethylene film)	
Pack factor	24 units per carton	
Packaging	Paper and film peel pouch 1 polyliner contains all the units in a carton box	
Manufacturer	Medline International France SAS EN ISO 13485	
Regulation	In compliance with Medical Devices Directive 93/42/EEC MDD Class: I sterile HP: High performance according to EN 13795 standard AAMI level : 4	



© 2015 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_Gown01_EN01

01.01.INT.12