



**PROXIMA FLUID COLLECTION POUCH
DYNJSD1016
TECHNICAL DATA**

Item Number	Dimensions (cm)	Qty / Case	Color	Sterilization
DYNJSD1016	48 x 50	40	Blue	Sterile

Description : PROXIMA FLUID COLLECTION POUCH

- ✓ Fluid Pouch 50x48 cm with adhesive, Non linting, tear resistant, poly release liner, moldable strip, micro screen for improved filtering system, suction port with cap, heavy low glare matte finish film.

Material : Polyethylene Film
Adhesive


Storage Conditions: Dry and Clean Conditions

Available Sterilization: Ethylene Oxide

Packaging: Single Sterile

Manufacturer: Conforms to ISO 13485

Law : Conforms to EU Directive 93/42/EEC
- Class I



Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.



Medline Industries, Inc
The Company of Choice
www.medline.com

 Medline
 One Medline Place
 Mundelein, Illinois 60060 USA
 Telephone +1 847 643 3292
 Fax +1 847 643 3295
International@medline.com