



**PROXIMA VERTICAL ISOLATION INCISE DRAPE**  
**DYNJSD1017**  
**TECHNICAL DATA**

Item Number	Dimensions (cm)	Qty / Case	Color	Sterilization
<b>DYNJSD1017</b>	213 x 317	15	Clear	Sterile

**Description :** PROXIMA VERTICAL ISOLATION INCISE DRAPE

- ✓ Tube Holders
- ✓ Instrument Pouches
- ✓ Incise Area 48x25cm
- ✓ Fluid Control Pouch

**Material :** SMS Polypropylene  
 Polyethylene Film  
 Incise


**Storage Conditions:** Dry and Clean Conditions

**Available Sterilization:** Ethylene Oxide

**Packaging:** Single Sterile

**Manufacturer:** Conforms to ISO 13485  
 Conforms to EN 13795 – High Performance

**Law :** Conforms to EU Directive 93/42/EEC  
 - Class I



**317 cm x 213 cm**  
**125 in. x 84 in.**

Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.



**Medline Industries, Inc**  
**The Company of Choice**  
[www.medline.com](http://www.medline.com)


**Medline**  
 One Medline Place  
 Mundelein, Illinois 60060 USA  
 Telephone +1 847 643 3292  
 Fax +1 847 643 3295  
[International@medline.com](mailto:International@medline.com)

DYNJSD1017 – RA013  
 Effective: 06/15/11

CO# C11-00915  
 Approved by: P. Nichol, D. Podlesak, P. Vasquez

Originator: N. Westburg  
 Released by: D. Moore