




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/en/](http://www.medline.com/en/)



## Technical sheet

# Cesarean Pouch

### Reference(s)

**DYNJSD2057**

### Description

360 degree cesarean pouch 91x76 cm, with incise film 38x36 cm, 2 cord holders

### Main materials

Matte finish polyethylene film

### Packaging

Paper and film peel pouch  
10 units per case

### Manufacturer

Medline Industries Inc Mundelein USA  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile



Latex free



0086



STERILE EO

© 2012 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-DYNJSD2057-  
EN01

01.01.INT.12