




**Medline Industries Ltd**  
 No 5 Booths Park - Chelford Road  
 Knutsford – Cheshire - WA16 8GS – England  
 Phone: 0844 334 5237  
 Fax : 0844 334 5238  
[www.medline.com/uk](http://www.medline.com/uk)



# Technical sheet

# PERFORATED DRAIN

## Medline Closed Wound Drainage

### References & Description


Medline Perforated Drains are used to remove exudate and excess fluid from a patient's surgical site wound to promote healing. The perforated end of the drain is placed in the surgical wound site to be enclosed within the patient. The distal end of the drain tubing is connected to a reservoir bulb for fluid collection.

#### Medline Perforated ROUND Drain:

- Are radiopaque stripe running the entire length of the drain tubing for x-ray detection.
- Come in a one-piece design and are hubless.
- Contain two white indicators points to aid in accurate drain placement.
  - o First (one point) : 5 cm after the perforated end of the drain
  - o Second (two points) : 10 cm after the perforated end of the drain

Total length: 125 cm

- A. Perforation area (22 cm)
- B. External transparent part (103 cm)
- C. Trocar is curved on the top part (15 cm)




Item codes	Diameter	Trocar	Connector / Adaptor
DYNJWE1320	7Fr	No	Yes (individual pouch)
DYNJWE1321A	10Fr		No
DYNJWE1323A	15Fr		No
DYNJWE1325A	19Fr		No
DYNJWE0321	10Fr	Yes (10 Fr)	Yes (individual pouch)
DYNJWE0323	15Fr	Yes (15 Fr)	No
DYNJWE0325	19Fr	Yes (19 Fr)	No

#### Medline Perforated FLAT Drain:

- Manufactured with barium sulfate for x-ray detection of drain
- have a two-piece design with inner lumen ribbing to help prevent drain collapse and clogging and help preserve drain patency

Total length: 101 cm

- A. Perforation area (3/4 = 14,5 or 4/4 = 19,5 cm)
- B. Internal white part (22 cm)
- C. External transparent part (79 cm)
- D. Trocar is curved on the top part (15 cm)



Item codes	Diameter	Perforation	Trocar	Reservoir
DYNJWE1308A	7mm	3/4	No	No
DYNJWE1309	10mm			
DYNJWE1310	7mm	4/4	Yes (15 Fr)	
DYNJWE1311	10mm		Yes (19 Fr)	
DYNJWE1410	7mm	3/4	No	Yes (100 mL)
DYNJWE1411	10mm			
DYNJWE1348	7mm	4/4	No	
DYNJWE1349	10mm			
DYNJWE1360	7mm	4/4	No	
DYNJWE1361	10mm			



© 1017 Medline Industries Inc. or one of its subsidiaries.  
 All rights reserved.

LIT-DrainPlaie\_Perf-EN02

01.01.INT.12

<b>Main materials</b>	<u>Drain</u> : silicone (Latex and DEHP free) - <u>Trocar</u> : Stainless steel
<b>Pack factor/Packaging</b>	Double sterile packaging (1 PE bag + 1 peel pouch) - 10 units / case Trocar is protected with a blue plastic protector.
<b>Sterilization</b>	Sterilized by Ethylene oxide
<b>Manufacturer</b>	DEGANIA SILICONE LTD - Degania Bet 1513000, Israel – EN ISO 13485 Distributed by Medline International France SAS
<b>Regulations</b>	In compliance with Medical Devices Directive 93/42/EEC MDD Class: IIb-Sterile



© 1017 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-DrainPlaie\_Perf-EN02

01.01.INT.12