




Medline International France SAS  
5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/en](http://www.medline.com/en)



# UltraGard™

# Surgical Gown

**Reference(s)**

Reference	Length (cm)	Size	Hand towels	Unit/case
E3505CEA	115	S/M	2	36
E3506CEA	120	L	2	36
E3516CEA	150	XLL	2	32
E3545CEA	130	XL	2	36
FE3545CE	140	XL	2	32
E3585CEA	160	2XLL	2	32

**Description**

A SMS wrap 60x60 cm contains:  
- 1 SMMMS sterile back gown, hook and loop adjustable neck - (EN13795-1 SP)  
- 2 absorbent hand towels, 30x39 cm unfolded

**Main materials**

SMMMS (polypropylene)

**Pack factor**

Paper and film peel pouch  
2 polyliners contain all the units in a carton box

**Manufacturer**

Medline International France SAS  
EN ISO 13485

**Regulation**

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile

- Please refer to complete labeling and instruction for use before using the product.



© 2020 Medline Industries Inc. or one of its subsidiaries. All rights reserved.  
TDS\_GownUltraGardStandard01\_EN05