




Medline International France SAS

5 rue Charles Lindbergh  
44110 Châteaubriant

[www.medline.com/en/](http://www.medline.com/en/)



## OPS™ Essential

## Arthroscopy pack

### Reference(s)

**ES10129CE**

### Description

A SMS wrap 90x90 cm contains:

- 1 back table cover, 140x190 cm, 50 µm containing :
  - 1 bi-laminate drape sheet, 150x150 cm
  - 1 bi-laminate arthroscopy T drape, 221/290x320 cm, 5 cm diameter elastic fenestration, armboard covers, fluid collection pouch, 2 hook & loop tube holders

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)

### Pack factor

10 units per carton

### Packaging

Vented sterilization pouch  
1 polyliner contains all the units in a carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I sterile  
HP: High performance according to EN 13795 standard



Latex free



STERILE EO



© 2014 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-ES10129CE-EN01

01.01.INT.12