




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Shoulder drape, beach chair

Reference(s)

ES10369CE

Description

A SMS wrap 90x90 cm contains:

- 1 bi-laminate shoulder drape, 411x262 cm, 10x10 cm pear shaped elastic fenestration, fluid collection pouch, 2 hook and loop tube holders

Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)

Pack factor

6 units per carton

Packaging

Vented sterilization pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795-1 standard



Latex free



2797



STERILE EO



© 2020 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_OrthopedicDrape_ES10369CE_EN01

01.INT.03.30