




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ ESSENTIAL

TUR Pack

Reference(s)

ES10381CE

Description

A sterile wrap 90x90 cm contains:

- 1 back table cover, 140x190 cm, 40 µm containing:

- 1 bi-laminate TUR drape, 250/280x280cm, fluid collection pouch with exit port, 6 cm diameter adhesive fenestration, 22x31 cm PE reinforcement with finger cot and 5 cm diameter fenestration
- 1 adhesive tape, 9x50 cm
- 2 absorbent towels, 30x20 cm unfolded

Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)

Pack factor

7 units per carton

Packaging

Vented sterilization pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



STERILE | EO



© 2013 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-ES10380CE-EN04

01.01.INT.12