




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Laparoscopy pack

Reference(s)

ES10496CE

Description

A SMS wrap 112x112 cm contains:

- 1 back table cover, 140x190 cm, 50 µm containing :
 - 1 bi-laminate Laparoscopy drape, 255x306 cm, 26x32 cm abdominal adhesive fenestration, 15x18 cm perineal adhesive fenestration with flap, leggings, 2 hook and loop tube holders
 - 1 clear plastic instrument pouch, 61x42,5 cm, 2 compartments, adhesive
 - 1 adhesive tape, 9x50 cm
 - 1 hook and loop tube holder 2.5x23 cm
 - 2 hand towels, 58x39 cm unfolded

Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)

Pack factor

7 units per carton

Packaging

Vented sterilization pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



Latex free



STERILE | EO

© 2015 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_GENSURGPACK_ES10496CE_EN01

01.01.INT.12