




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Drape sheet

Reference(s)

ES15226CE

Description

A peel pouch, 15x28 cm, contains:
- 1 bi-laminate drape 75x45 cm

Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)

Pack factor

320 units per carton
80 units per dispenser box

Packaging

Paper and film peel pouch
4 dispenser boxes contain all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



Latex free



0086



STERILE EO



© 2014 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-ES15226CE-EN03

01.01.INT.12