




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Universal pack

Reference(s)

ES29105CE

Description

A SMS wrap 90x90 cm contains:

- 1 back table cover, 140x190 cm, 50 µm containing:
 - 4 absorbent towels, 30x20 cm unfolded
 - 1 hook and loop tube holder, 14x2.5 cm
 - 1 adhesive tape, 9x50 cm
 - 2 bi-laminate side drapes 110x95 cm, 67x30 cm reinforcement, adhesive
 - 1 bi-laminate foot drape 200x200 cm, 67x30 cm reinforcement, adhesive
 - 1 bi-laminate head drape 260x170cm, 67x25 cm reinforcement, adhesive

Main materials

Absorbent and impervious bi-laminated material (polyethylene + polypropylene)
Absorbent Spunbond polypropylene reinforcement

Pack factor

8 units per carton

Packaging

Vented sterilization pouch
1 polyliner contains all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
HP: High performance according to EN 13795 standard



Latex free



0086



STERILE EO



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_UnivPack_ES29105CE_EN04

01.INT.03.30