




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



OPS™ Essential

Drape Sheet - Extra Large

Reference(s)	ES9444CE
Description	A SMS wrap 76x76 cm, contains: - 1 bi-laminate drape, 196x249 cm
Main materials	Absorbent and impervious bi-laminated material (polyethylene + polypropylene)
Pack factor	18 units per carton
Packaging	Vented sterilization pouch 1 polyliner contains all the units in a carton box
Manufacturer	Medline International France SAS EN ISO 13485
Regulation	In compliance with Medical Devices Directive 93/42/EEC MDD Class: I Sterile HP: High performance according to EN 13795 standard



© 2017 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_GenSurgDrape_ES9444CE_EN01

01.01.INT.12