



Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant
www.medline.com/en/



Technical sheet

Electrosurgical Pencil

Reference(s)

ESPB4001LQ

Description

Electrosurgical pencil with push button, stainless steel blade
Shaft diameter: 2,4 mm

Cable length Connector

5 m
International

Maximum electrical capacity

5 kVp

Main materials

Electrode/Blade: Stainless steel
Pencil: ABS 100

Pack factor

80 units per carton / 2 boxes per carton

Packaging

PE bag
Peel pouch

Manufacturer

Dongguan QueenMed Equipment Mfg. Co., Ltd. - China
EN ISO 13485:2012

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: IIb sterile
Electrical Safety Standard : IEC 60601-2-2



STERILE R

© 2017 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-ESPB4001LQ-EN01

01.01.INT.12