



Medline International France SAS
 5 rue Charles Lindbergh
 44110 Châteaubriant
 www.medline.com/en



Technical sheet

Skin Stapler

Reference(s)

F-35WB – wide staples
F-35RB – regular staples

Description

Disposable Skin Stapler with 35 stainless steel staples

Dimensions closed staples	Width	Height
Wide staples	7.2 mm	4.9 mm
Regular staples	5.9 mm	3.9 mm

Main materials

Body: Acrylonitrile butadiene styrene (ABS)
 Staples cartridge: Polycarbonate
 Staples: Medical-grade stainless steel

Pack factor

5 units per dispenser
 100 units per carton (20 dispenser boxes)

Packaging

Peel packaging

Manufacturer

Ningbo Advan Electrical Co.,Ltd
 EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
 MDD Class: IIa sterile



STERILE EO

© 2015 Medline Industries Inc. or one of its subsidiaries.
 All rights reserved.

LIT-SkinStapler-
 EN03

01.01.INT.12