




Medline Industries Ltd
3rd Floor, Quayside Wilderspool Business Park
Greenalls Avenue
Warrington WA4 6HL

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical Sheet

Foley Anchor

Reference(s)

Reference	Width (cm)	Length (cm)	Recommended usage	Pictures
ATA24ZT	10,5	6,1	Urinary Catheter Drainage Line	
FCS200ZT	11,4	7,6		

Description

Foley Anchor is a urinary catheter securement device that removes the frustration factor often experienced by caregivers and patients attempting to keep a Foley catheter securely in place.

Foley Anchor is easy to apply and remove and fits virtually any size catheter. The long wear time reduces nursing time and it prevents catheter dislodgement.

Main materials

Liner: Silicone treated paper
Securement band :Silicone
Medical grade adhesive
Base : PVC (FCS200ZT) / PE (ATA24ZT) foam

Pack factor

50 units per carton

Packaging

Pouch - A carton box contains all the units

Manufacturer

Centurion Medical Products Corporation (a Medline company) - EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile

- Please refer to complete labeling and instruction for use before using the product.

Shelf Life

5 years



© 2020 Medline Industries Inc. or one of its subsidiaries. All rights reserved.
TDS_Foley Anchor_EN02