



Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



Technical sheet

Clip Cap

Reference(s)

Reference	Color
FS65450	Blue
FS65460	White

Description

Clip cap, one size fits most

Main materials

Spunbond polypropylene, 14 gsm

Pack factor

1200 units per case

Packaging

12 polyliners contain all the units in a carton box

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC
MDD Class: 1 Non-Sterile



Latex free



© 2017 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_Headwear_FS654xx_EN02

01.01.INT.12