




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



Technical sheet

Anti-Skid Instrument Mat

Reference(s)

GSASM1CE

Description

Gold Standard Anti-Skid Instrument Mat, 25 x 41 cm
2 adhesive strips on the back for fixation

Main materials

Polyvinyl chloride foam

Pack factor

80 units per carton

Packaging

Peel pouch

Manufacturer

Medline Industries, Inc. USA
EN ISO 13485

Regulation

MDD Class: n/a



Latex free



STERILE | EO

© 2012 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LITAntiSkidInstrumentMat-EN01

01.01.INT.12