




Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en/



Technical Sheet

Sponge bag

Reference(s) IBE5242, IBE5637, IBE5652

Description

Reference	Dimensions (cm)	Number of compartments
IBE5242	52x42	40
IBE5637	56x37	6
IBE5652	56x52	20

Main materials Clear film (Polyethylene)

Pack factor 50 units per carton

Packaging Paper and film peel pouch
1 polyliner contains all the units in a carton box

Manufacturer Medline International France SAS
EN ISO 13485

Regulation In compliance with Medical Devices Directive 93/42/EEC
MDD Class: I sterile
In compliance with the standard EN13795



Latex free



0086



STERILE EO

© 2018 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS_SpongeBag_EN01

01.01.INT.12