



Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant, France

www.medline.com/en/

Technical sheet

Sterile White Lap Sponges 5's

Reference(s)

Item Codes	Size	Ply	Packaging	Pack Factor
LW5P203004	20cm x 30cm	4	5 pc/pk	30 pk/cs
LW5P454504	45cm x 45cm	4	5 pc/pk	20 pk/cs
LW5P454506	45cm x 45cm	6	5 pc/pk	20 pk/cs

Description

- X-ray detectable lap sponges, 100% cotton, prewashed, white, soft and absorbent, with a loop sewed into the sponge. Single use. Sterile. Do not use if the packaging (peel pouch) is opened or damaged. Do not resterilize. Two separate adhesive counting cards included.
- **Intended use:** For absorption of body fluids in or around the body and to retain organs during surgery

Main materials

Cotton – Type 20 Thread
Barium sulfate strands

Pack factor

Refers to "References" item

Packaging

Double Soft Pack (paper-plastic)
Carton box

Manufacturer

Medline Industries, Inc.
One Medline Place
Mundelein, Illinois 60060 (USA)
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: IIa Sterile

In compliance with EN 14079



STERILE EO



© 2013 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

LIT-TWC-LW5P-EN01

01.01.INT.12