



Medline Industries Ltd  
3rd Floor, Quayside,  
Wilderspool Business Park  
Greenalls Avenue  
Warrington WA4 6HL

Phone : +44 (0)800 151 3587

[www.medline.eu/en](http://www.medline.eu/en)



## Technical Sheet

# Bordered Gauze Dressing

### Reference(s)

Item Code	Dressing Size (cm)	Absorbent Pad Size (cm)
MSCE3223	5 x 7	3 x 4
MSCE3278	6.5 x 8	3.5 x 5
MSCE3246	10 x 15	5 x 10
MSCE32410	10 x 25	5 x 20
MSCE32414	10 x 35	5 x 30
MSCE3244	10 x 10	6 x 6
MSCE3248	10 x 20	5 x 15

### Description

Bordered gauze dressings come in different sizes to meet the specific wound or incision size.

- Deluxe soaker pad
- Non-woven adhesive border
- Water resistant backing
- Ideal for incision sites

### Main materials

Polyester and viscose

### Pack factor

15 dressings per box, 10 boxes per case (150 dressings per case)

### Manufacturer

Changzhou Hualian Health Dressing Co. - EN ISO 13485  
EC Representative: Shanghai International Holding Corp. GmbH, Eiffestrasse 80,20537  
Hambourg, Germany  
Distributed by Medline

### Regulation

In compliance with Medical Devices Directive 93/42/EEC

MDD Class: I Sterile

- Please refer to complete labeling and instruction for use before using the product.



© 2020 Medline Industries Inc. or one of its subsidiaries. All rights reserved.  
TDS\_Bordered-Gauze\_EN02