



## ULTRASORBS AP DISPOSABLE UNDERPAD TECHNICAL DATA

| AIR PERMEABLE |              |                     |                    |                       |            |                    |
|---------------|--------------|---------------------|--------------------|-----------------------|------------|--------------------|
| Product No.   | Size         | Backing Print Color | Packaging          | Case Quantity         | Absorbency | Box Dimensions     |
| ULTRASRBE6191 | 61 x 91 cm   | Blue                | Poly Bag/Corrugate | 70; 7 bags of 10      | 1800-2300g | 17.75x8.5x17.50 in |
| ULTRASRBE2540 | 25 x 40 cm   | Blue                | Poly Bag/Corrugate | 100; 10 bags of 10    | 200-300g   | 17.3x10.6x5.5in    |
| ULTRASRBE4561 | 45 x 61 cm   | Blue                | Poly Bag/Corrugate | 60; 6 bags of 10      | 900-1150g  | 17.7x9.8x7.9in     |
| USAPE6191UF   | 61 x 91 cm   | Blue                | Poly Bag/Corrugate | 70; 7 bags of 10      | 1800-2300g | 17.75x8.5x17.50 in |
| USAPE6191ES   | 61 x 91 cm   | Green               | Poly Bag/Corrugate | 60; 6 bags of 10      | 1800-2300g | 17.75x8.5x17.50 in |
| MSCE4090PUPS  | 102 x 229 cm | Blue                | Poly Bag/Corrugate | 30; 1 bag of 30       | 4000-4500g | 24x16x13 in        |
| MSCE4090PUPSB | 102 x 229 cm | Blue                | Poly Bag/Corrugate | 30; 30 bags of 1 each | 4000-4500g | 24x16x13 in        |

| NON-AIR PERMEABLE |            |     |                    |                  |            |           |
|-------------------|------------|-----|--------------------|------------------|------------|-----------|
| MSCE282100        | 61 x 91 cm | N/A | Poly Bag/Corrugate | 70; 7 bags of 10 | 1800-2300g | 18x9x18in |

- Description :** A package contains : Disposable Underpad
- Material :** Spun-lace polypropylene, Airlaid core with Super absorbent polymer, Air permeable BTBS film/nonwoven Backsheet (on all except MSCE282100)
- Storage Conditions:** General Warehouse Conditions, Keep Dry
- Sterilization :** NON-STERILE
- Manufacturer :** Conforms to EN ISO13745
- Law :** Conforms to EU Directive 93/42/EEC- Class I



**Medline Industries, Inc**  
The Company of Choice  
[www.medline.com](http://www.medline.com)

 Medline  
One MedlinePlace  
Mundelein, Illinois 60060 USA  
Telephone +1 847 643 3292  
Fax +1 847 643 3295  
[International@medline.com](mailto:International@medline.com)

Ultrisorb – RA013  
Effective: 07/24/12

CO# C12-01631  
Approvers : H. Hood, P. Vasquez, S. Smith

Originator : S. Smith  
Released by : W. O'Connor