

# PROHIBIT® SERIES FACE MASKS

## Procedural Mask – Product Data Sheet

**Reference:** NON27378/NONE27378

**Description:** Procedural Face Mask with Earloops, Blue



**Materials:**

<b>Inner facing:</b>	Polypropylene, White
<b>Outer Facing:</b>	Polypropylene, Blue
<b>Filter Media:</b>	Meltblown Polypropylene, White
<b>Loops:</b>	Thermal Bonded Polypropylene, White

**Packaging:** Packaged in Bulk, 50 Masks/Carton, 6 Cartons/Case

**Manufacturer:** Conforms to EN ISO 13485, Conforms to EN14683, Type II

**Law:** Conforms to EU Directive 93/42/EEC - Class I non-Sterile

**Environment:** Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.

Physical Properties	EN14683 Type II Requirements	EN14683 Medline Results
Bacterial Filtration Efficiency (BFE%) at 3 micron	≥ 98%	99.7%
Particle Filtration Efficiency (PFE%) 0.1 micron	n/a	97.6%
Differential Pressure (Delta P)	< 3.0 mmH2O/cm2	2.52 mmH2O/cm2

*Tests have been performed by: Nelson Laboratories, 6280 South Redwood Road, Salt Lake City, UT 84123*

©2011 Medline Industries, Inc., Mundelein, Illinois 60060 USA.  
Medline is a registered trademark of Medline Industries, Inc



Medline  
One Medline Place  
Mundelein, Illinois 60060 USA  
Telephone +1 847 643 3292  
Fax +1 847 643 3295  
[International@medline.com](mailto:International@medline.com)



[www.medline.com](http://www.medline.com)