



Medline International France SAS  
 5 rue Charles Lindbergh  
 44110 Châteaubriant  
[www.medline.com/en/](http://www.medline.com/en/)



## Technical Sheet

## Procedure Facemask

### Reference(s)

**NONE27378**

### Description

Procedure facemask with earloops, blue

### Main materials

Outer layer: Spunbond (Polypropylene)  
 Filter layer: Meltblown (Polypropylene)  
 Inner layer: Spunbond (Polypropylene and Polyethylene)  
 Earloops: Latex-free elastic

### Pack factor

300 units per carton

### Packaging

6 dispensers contain all the units in a carton box

### Manufacturer

Medline International France SAS  
 EN ISO 13485

### Regulation

In compliance with Medical Devices Regulation 2017/745 - MD Class: I Non Sterile  
 In compliance with EN14683 – Type II

### Performance

Physical Properties	EN14683 Type II Requirements	EN14683 Medline Results
Bacterial Filtration Efficiency (BFE%) at 3 micron	≥ 98	99.6
Differential Pressure (Pa/cm <sup>2</sup> )	< 29.4	25.3
Microbial cleanliness (cfu/g)	≤30	5



© 2020 Medline Industries Inc. or one of its subsidiaries.  
 All rights reserved.

TDS\_FaceMask\_NONE27378\_EN02 01.INT.03.30