

PROHIBIT® SERIES FACE MASKS

Procedural Mask – Product Data Sheet

Reference: NON27386/NONE27386
Description: Surgical Face Mask, Hypoallergenic



Materials:
Inner facing: Thermal Bond Polypropylene, White
Outer Facing: Polypropylene, White
Filter Media: Meltblown Polypropylene, White
Ties: Thermal Bonded Polypropylene, White

Packaging: Packaged in Bulk, 50/box, 6 box/cs, (300ea/cs)

Manufacturer: Conforms to EN ISO 13485, Conforms to EN14683, Type II

Law: Conforms to EU Directive 93/42/EEC - Class I non-Sterile

Environment: Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.

Physical Properties	EN14683 Type II Requirements	EN14683 Medline Results
Bacterial Filtration Efficiency (BFE%) at 3 micron	≥ 98%	99.8%
Particle Filtration Efficiency (PFE%) 0.1 micron	n/a	n/a
Differential Pressure (Delta P)	< 3.0 mmH2O/cm2	2.86 mmH2O/cm2

Tests have been performed by: Nelson Laboratories, 6280 South Redwood Road, Salt Lake City, UT 84123

©2011 Medline Industries, Inc., Mundelein, Illinois 60060 USA.
Medline is a registered trademark of Medline Industries, Inc

 Medline
One Medline Place
Mundelein, Illinois 60060 USA
Telephone +1 847 643 3292
Fax +1 847 643 3295
International@medline.com



www.medline.com

NON27386, NONE27386 – RA013
Effective 08/05/11

CO# C11- 01233
Approved by: B. Tompkins, P. Vasquez, P. Nichol

Originator: N. Westburg
Released by: W. O'Connor