



NON-STERILE

ISOLATION GOWNS
TECHNICAL DATA

Product No.	Size	Weight	Length	Packaging	Quantity	Case Dimensions
CRIE4000	Regular/ Large	20 grams SPP	45 inches (minimum)	Inner Poly Bag Corrugate Case	50/cs	14.7Lx11.9Wx10.2 H
CRIE4001	X-Large	20 grams SPP	50 inches (minimum)	Inner Poly Bag Corrugate Case	50/cs	15.0Lx12.75Wx10. 75H
CRIE4000G	Regular/ Large	20 grams SPP	45 inches (minimum)	Inner Poly Bag Corrugate Case	50/cs	14.1Lx11.9Wx9.45 H
CRIE4010CF	Regular/ Large	20 grams SPP	45 inches (minimum)	Inner Poly Bag Corrugate Case	50/cs	14.75Lx12.0Wx13. 0H
CRIE4011CF	X-Large	20 grams SPP	50 inches (minimum)	Inner Poly Bag Corrugate Case	50/cs	15.0Lx11.8Wx13.0 H
NONE27SMS3	Regular/ Large	25 grams SMS	45 inches (minimum)	Inner Poly Bag Corrugate Case	100/cs	14.9Lx12.0Wx16.6 H
NONE27SMS3XL	X-Large	25 grams SMS	50 inches (minimum)	Inner Poly Bag Corrugate Case	100/cs	15.35Lx13.0Wx18. 5H
CRIE5000	Regular/ Large	25 grams PE	45 inches (minimum)	Individual Poly Bag, Inner Box, Corrugate Case	15/bx, 75/cs	Box 6.5Lx4.0Wx8.25H Case 21.5Lx8.25Wx9.25H
CRIE5001	X-Large	25 grams PE	53 inches (minimum)	Individual Poly Bag, Inner Box, Corrugate Case	15/bx, 75/cs	Box 6.5Lx4.0Wx9.0H Case 21.5Lx7.75Wx10.25H
NONE27116	Regular/ Large	48 grams Coated SPP PE	45 inches (minimum)	Inner Poly Bag, Corrugate Case	50/cs	Case 14.56Lx11.81Wx12.60H
NONE27116XL	X-Large	48 grams Coated SPP PE	50 inches (minimum)	Inner Poly Bag, Corrugate Case	50/cs	Case 14.56Lx11.81Wx12.60H
NONE27457	Regular/ Large	82 grams Microporous Rubber B- Component Polypropylene Material	45 inches (minimum)	Inner Poly Bag, Corrugate Case	30/cs	Case 15Lx12.5Wx13H
NONE27457XL	X-Large	82 grams Microporous Rubber B- Component Polypropylene Material	50 inches (minimum)	Inner Poly Bag, Corrugate Case	30/cs	Case 15.8Lx13.5Wx11.4H
NONE27457XXL	XX-Large	82 grams Microporous Rubber B- Component Polypropylene Material	54 inches (minimum)	Inner Poly Bag, Corrugate Case	30/cs	Case 17.2Lx12.5Wx14.9H

Description : Gown protects healthcare worker and patient

Material : Gown made from
100% Spunbond Polypropylene Nonwoven Material
100% SMS Polypropylene Nonwoven Material
Coated Spunbond Polypropylene Polyethylene or
Polyethylene (PE) Material

Storage Conditions: At room temperature

Packaging: Ten gowns per poly bag (unsealed), 5 or 10 bags per case or
15 gowns per inner box, 5 boxes per case

Sterilization: NON-STERILE

Manufacturer : Conforms to ISO 13485

Law : Conforms to EU Directive 93/42/EEC
- Class I



Medline Industries, Inc. is committed to reducing environmental impact through increased use of post-consumer recycled materials.

Medline Industries, Inc
The Company of Choice
www.medline.com

 Medline
One Medline Place
Mundelein, Illinois 60060 USA
Telephone +1 800 MEDLINE
International@medline.com

NON CRI Gowns – RA012
Effective: 05/14/12

CO# C12-00972
Approvers : S. Pasko, N. O'Connor, P. Vasquez

Originator : 05/14/12
Released by : W. O'Connor