



Medline International France SAS

5 rue Charles Lindbergh
44110 Châteaubriant

www.medline.com/en



Technical sheet

Footwear

Reference(s)

Reference	Description	Color	Size	Main materials	Basis weight	Packaging
CRIE2002	Shoe covers, non-skid bottom	Blue	L	SMS (polypropylene)	35g/m ²	100/bx 300/cs
CRIE2010	Shoe cover, plain bottom	Blue	One Size Fits Most	Polyethylene	27g/m ²	1000/cs
NONE28858	Shoe cover, non skid bottom	Blue	L	SMS (polypropylene)	35g/m ²	100/bx, 300/cs
NONE28859	Shoe cover, non skid bottom	Blue	XL	SMS (polypropylene)	35g/m ²	100/bx, 200/cs
NONE27143XL	Boot cover, non skid bottom, elastic band at the ankle	Blue	XL	Polyethylene coated on spunbond polypropylene	50g/m ²	50/bx, 150/cs
NONE27144	Boot cover, non-skid pattern, hook and loop strap	Blue	L	Polyethylene coated on spunbond polypropylene	50g/m ²	50/bx, 150/cs
NONE27144XL	Boot cover, non-skid pattern, hook and loop strap	Blue	XL	Polyethylene coated on spunbond polypropylene	50g/m ²	50/bx, 150/cs
NONE27348PXL	Boot cover, foam non-skid strips, hook and loop strap	Blue	XL	Spunbond polypropylene and polyethylene	82g/m ²	50/bx, 150/cs

Manufacturer

Medline International France SAS
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/EEC and Personal Protective Equipment 89/686/EEC

MDD Class: 1 Non-Sterile

PPE Class: 1



© 2013 Medline Industries Inc. or one of its subsidiaries. All rights reserved.

LIT-Footwear-EN02

01.01.INT.12