



Medline International France SAS
5 rue Charles Lindbergh
44110 Châteaubriant, France

www.medline.com/en



Technical Sheet

Sterile Aspiration Tubing

Reference(s)

Reference	Description	Qty/case
ORNEX56A	Tubing, sterile, 1.80 m – internal diameter 5 mm	50
ORNEX66A	Tubing, sterile, 1.80 m – internal diameter 6 mm	45
ORNEX76A	Tubing, sterile, 1.80 m – internal diameter 7 mm	40
ORNEX510	Tubing, sterile, 3.0 m – internal diameter 5 mm	30
ORNEX610	Tubing, sterile, 3.0 m – internal diameter 6 mm	25
ORNEX710	Tubing, sterile, 3.0 m – internal diameter 7 mm	20
ORNEX712	Tubing, sterile, 3.60 m – internal diameter 7 mm	20

Description

- Sterile aspiration tubing – Patient/canister.
Sterile aspiration tubing with universal rubber connectors on each end for a flawless patient canister connection and removable male connector.
Diameters available: 5 mm - 6 mm - 7 mm.
Tubing is rolled and packaged in individual double packaging blister.

Main materials

Transparent PVC, containing DEHP

Pack factor

1 each per pouch
Cf item " Reference(s)" for pack factor per case

Packaging

Sterilization pouch (polybag and paper/poly pouch)
carton box

Manufacturer

Medline Industries, Inc
EN ISO 13485

Regulation

In compliance with Medical Devices Directive 93/42/CEE
MDD Class : IIa Sterile



STERILE EO



© 2013 Medline Industries Inc. ou une de ses filiales.
Tous droits réservés.

LIT- ORNEXSterile-EN01

01.01.INT.12