




Medline Industries Ltd
 No 5 Booths Park
 Chelford Road
 Knutsford - Cheshire
 WA16 8GS – England
 Phone: 0844 334 5237
 Fax : 0844 334 5238

www.medline.com/en/



Technical sheet

YANKAUER SET Tubing + Yankauer

Reference(s)

ORNEX66AY1200, ORNEX610Y1800, ORNEX610Y2200, ORNEX610Y80

Description

Sterile set of suction tubing with 2 female ribbed connectors and Yankauer pre-attached.

Code	Description	Pack factor (each/case)
ORNEX66AY1200	Suction set Yankauer with rigid cannula with lateral holes (2) and beveled tips for controlled atraumatic suction. Yankauer length (in mm) : 275mm (±15mm) Yankauer diameter (in CH): 12CH Tubing length (in m) : 1.80m Tubing inner diameter (in mm) : 6mm	30
ORNEX610Y1800	Suction set Yankauer with rigid cannula with lateral holes (6) and beveled tips for controlled atraumatic suction. Yankauer length (in mm) : 250mm (±15mm) Yankauer diameter (in CH): 18CH Tubing length (in m) : 3.00m Tubing inner diameter (in mm) : 6mm	30
ORNEX610Y2200	Suction set Yankauer with rigid cannula with lateral holes (6) and beveled tips for controlled atraumatic suction. Yankauer length (in mm) : 275mm (±15mm) Yankauer diameter (in CH): 22CH Tubing length (in m) : 3.00m Tubing inner diameter (in mm) : 6mm	30
ORNEX610Y80	Suction set Yankauer with rigid ergonomic cannula with lateral holes (4) and beveled tips for controlled atraumatic suction. Non-slip grip handle design for comfort and security during long procedures. Yankauer length (in mm) : 288mm (±20mm) Yankauer diameter (in CH): 12CH Tubing length (in m) : 3.00m Tubing inner diameter (in mm) : 6mm	30



© 2019 Medline Industries Inc. or one of its subsidiaries.
 All rights reserved.

TDS- YankauerOrnex-EN01

01.INT.03.30

Main materials	<u>Tubing</u> : Light blue PVC. Latex Free - <u>Connectors</u> : Blue PVC. Latex free Yankauer Y1200, Y1800, Y100 : PVC. Latex free. DEHP free Yankauer Y80: K-resin. Latex free. DEHP free
Pack factor / Packaging	Sterile double wrapped transparent individual poly bag and peel pouch – Pack factor: see above
Sterilization	Sterilized by Ethylene Oxide
Manufacturer	Gcmedica Enterprise Ltd. (Wuxi) – ISO 13485 Distributed by Medline International France SAS.
Regulations	In compliance with Medical Devices Directive 93/42/EEC - MDD Class: IIa sterile



© 2019 Medline Industries Inc. or one of its subsidiaries.
All rights reserved.

TDS- YankauerOrnex-EN01

01.INT.03.30